REF 899021 | 2022-03-24 © Sutter Medizintechnik GmbH

Lagerungstrays / Storage Trays

REF:

701724, 701725, 701740, 701746, 701747,
701764, 701765, 701766, 701767, 701775-01,
701775-02, 701775-05, 701775-10, 701778-01,
701778-02, 701778-05, 701778-10

DE Lagerungstrays (Lagerungstrays aus Edelstahl)	5	FI Varastointilokero (Varastointilokerot ruostumattomasta teräksestä)	31
EN Storage Trays (Stainless steel storage trays)	8	NO Lagringsbrett i rustfritt stål	33
FR Plateau de Rangement (Plateau de rangement en acier inoxydable)	11	PL Tacki Magazynujące (Tacki magazynujące ze stali nierdzewnej)	36
ES Bandejas de Almacenamiento (Bandejas de almacenamiento de acero inoxidable)	14	BG Тави за Съхранение (Тави за съхранение от висококачествена стомана)	39
IT Vassoi Portastrumenti (Vassoi portastrumenti)	17	RO Tăvițe Depozitare (Tăvițe depozitare din inox)	42
NL Bewaarschotels (Bewaarschotels uit roestvrij staal)	20	SK Skladovacie Podnosy (Skladovacie podnosy z ušľachtilej ocele)	45
PT Bandejas de Colocação (Bandejas de colocação de aço inoxidável)	23	SL Pladnji za Shranjevanje (Pladnji za skladiščenje iz nerjavnega jekla)	47
SV Lagringstråg (Lagringstråg i rostfritt)	26	CZ Úložné Boxy (Úložné boxy z nerezové oceli)	50
DK Opbevaringsbakker (Opbevaringsbakker af rustfrit stål)	28	HU Tárolótálcák (Nemesacél tárolótálcák)	52

LT	Laikymo Padėklai (Laikymo padėklai iš nerūdijančiojo plieno)	55
LV	Glabāšanas Paplātes (Glabāšanas paplātes no nerūsējošā tērauda)	57
EL	Δίσκοι Αποθήκευσης (Δίσκοι αποθήκευσης από ανοξείδωτο ασάλι)	60
HR	Pladnjevi za spremanje od nehrđajućeg čelika	63
RU	Лотки (Лотки из нержавеющей стали)	65
SR	Kasete od nerđajućeg čelika	68
MK	Фиоки за складирање не’рѓосувачки челик	71
TR	Depolama Tekneleri (Paslanmaz çelikten oluşan depolama tekneleri)	74

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STERRAD® is a trademark of Johnson & Johnson, Inc.

Product / User / Disposal:

Medical devices may only be used and disposed of by qualified medical personnel!
Notwithstanding these instructions, reading the instructions for use provided with the stored instruments and other accessories being used is mandatory.

Intended purpose:

Sutter storage trays are intended for the safe storage of suitable Sutter instruments during reconditioning in a cleaning and disinfection device (CDD) and during sterilization in hot steam.

Service life:

A minimum of 100 reconditioning cycles can be expected with proper use.

Note:
For pre-cleaning, observe the instructions for use of the instruments being stored.
⚠ The instructions for use of the instruments being stored always take precedence over these instructions for use.

Machine cleaning / Disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

The fundamental suitability of the storage trays for cleaning instruments in a CDD was proven with the following program:

Program steps	Parameters
Pre-rinse with cold tap water	3 minutes
Clean with 0.5 % deconex® 28 ALKA ONE-x (Borer Chemie) at 70 °C	5 minutes
Mid-cycle and/or final rinse with warm tap water (40-45 °C)	1 minute
Rinse with deionized water	1 minute

- Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the trays and are permitted according to the A₀ concept, for instance thermal disinfection at 90 °C, 5 minutes, comparable A₀ value >3000. If a different cleaning agent is used,

only use products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie). Please contact your supplier or hygiene officer in case of questions.

⚠ Caution, instruments with internal lumina (such as irrigation forceps):

- Note special pre-cleaning requirements!
- Use connection options for flushing if possible!
- Please contact your hygiene officer in case of questions!

Inspection:

Prior to subsequent sterilization, perform a visual inspection and check for correct positioning, cleanliness, and integrity of the storage tray as well as the instruments it contains.

Packaging:

⚠ Package cleaned and disinfected storage trays in disposable sterilization packaging (single or double packaging), or place the storage tray with the cleaned and disinfected instruments into a suitable sterilization container meeting the following requirements:

- EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
- Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

Sterilization:

Storage trays / instruments must be cleaned and disinfected prior to sterilization.

The fundamental suitability of the storage trays for the sterilization of instruments in a steam sterilizer (according to EN 13060 and/or EN 285 and validated according to EN ISO 17665) was proven with the following program:

Program steps	Parameters
Process	3x fractionated vacuum process
Sterilization temperature	Min. 132 °C
Sterilization time (holding time at sterilization temperature)	3 minutes
Drying time	Min. 1 minute

⚠ Do not sterilize in hot air!

⚠ Do not sterilize in STERRAD®!

⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the instruments and storage tray and do not use it again.

Maintenance:

In general, do not use metal brushes, steel wool, scouring sponges, or abrasive cleaning agents to clean Sutter instruments, Sutter sterilization trays, and Sutter sterilization containers.

Do not expose any of the products to temperatures in excess of 141 °C.

Storage / Transportation:

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging.

For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Please note:

The instructions listed above have been validated by the manufacturer as suitable for preparing a medical device for reuse. It is the responsibility of the reconditioner to ensure that actual reconditioning with the equipment, materials, and personnel used in the reconditioning facility achieves the desired result.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The product may only be repaired by the manufacturer or an agent expressly authorized by the manufacturer. Otherwise the warranty and any possible additional liability claims against the manufacturer are voided.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik GmbH.

Changes reserved. Current version available at www.sutter-med.de.

STERRAD® is a trademark of Johnson & Johnson, Inc.

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