TAB 1			REF
A		-> Sutter BM-780 II, Martin, Berchtold, Aesculap, Erbet-Serie	36 01 85
B		-> Erbe ACC / IC / ViO, Karl Storz	36 01 86
C		-> Sutter CURIS®, Valleylab, Conmed, Bowie, Bowa	36 01 87
D		992901032	71 50 10
E		992901022	71 50 15
F		993801022	71 50 17
G		992901032	71 50 18
H		992901032	71 50 19

Deutsch

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STERRAD® is a trademark of Johnson & Johnson, Inc.

Product / User / Disposal:
Electrosurgery accessories may only be used and disposed of by qualified medical personnel!
Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.

⚠ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

Intended use:

Monopolar suction tube for suctioning smoke and fluids and for coagulating soft tissue during surgical procedures.

Service life:

A minimum of 20 reconditioning cycles can be expected with proper use.

Before using:

⚠ Inspect the product for cleanliness, intact insulation, and integrity before each use.
We recommend checking the insulation with a suitable test device.
⚠ Only use sterilized products that are in flawless condition!
A certain discoloration of the tips on non-stick instruments is normal and harmless.
Do not use in case of discernible damage to the electrical insulation, visible corrosion, or if the product exhibits nicks, dents, chipping, or cracks. Properly dispose of the product in such cases.
⚠ Only use the instrument with a suitable connection cable from Sutter! (see TAB1-A-C).
To connect the product, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!
⚠ Never use monopolar instruments without a neutral electrode! Verify that the neutral electrode is applied correctly.
⚠ Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory. Carefully follow the instructions provided by the manufacturer of the electrosurgery device being used, and read all warnings and safety notices.

During use:

Always work with the lowest output setting for the desired surgical effect.
Operation only in COAG mode with max. 40 watts. Brief, pulse-like activation if possible. Avoid continuous activation > 15 seconds.
Do not use in spray / fulguration mode!
Avoid spark formation!
⚠ After application, the tips may be so hot they cause burns!
⚠ Never lay down instruments on or in the immediate vicinity of the patient!
⚠ Position cables so the possibility of contact with the patient and other electrical lines is excluded.
⚠ When instruments are temporarily unused, store them so they are isolated from the patient.
⚠ Not for use with or in the vicinity of combustible or explosive substances!
⚠ When flushing instruments are used simultaneously, use non-conductive flushing fluid if possible!
⚠ Maximum allowable voltage 3000 Vp.

Reconditioning:

General information:

Observe national guidelines and regulations!
Disconnect the instrument from the cable!
Do not allow blood and tissue residues to dry up!
Remove blood and tissue residues with a soft cloth or brush!
Do not use sharp tools / scouring agents!
⚠ Do not immerse in hydrogen peroxide (H₂O₂)!

Manual cleaning and disinfection:

⚠ The instrument must always be reconditioned by machine – no manual cleaning!
According to the recommendations of the DGSV (“Deutsche Gesellschaft für Sterilgutversorgung”, German Association for Sterile Goods Supply), monopolar suction tubes are assigned to risk group B*. Machine cleaning is required for such products.

**This classification was assigned according to the DGSV 2013 flow chart for the classification of medical devices based on the recommendations of the KRINKO (“Krankenhausthygiene und Infektionsprävention Kommission”, Hospital Hygiene and Infection Prevention Commission)/BIARM (“Bundesinstitut für Arzneimittel und Medizinprodukte”, Federal Institute for Drugs and Medical Devices), Bundesgesundheitsblatt (Federal Health Bulletin) 2012; 55:1244-1310*

Manual pre-cleaning:

To make subsequent cleaning easier, do not allow tissue residue and blood to dry on during use.
Clean the instruments to remove major contamination directly after use (within 1 hour at most). Use running water or a disinfectant solution. The disinfectant must be aldehyde-free (to prevent fixing blood particles) with proven effectiveness (for instance VAH (“Verbund für Angewandte Hygiene e.V.”, Applied Hygiene Association)/DGHM (“Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.”, German Association for Hygiene and Microbiology) approval, CE marking) and be suitable for disinfecting the instruments (see “Material resistance”).

Only use a soft brush, soft cloth to manually clean or manually clean – no metal brushes, steel wool, scouring sponges, or abrasive cleaning agents.

Flush the suction tube five times with a disposable syringe (minimum volume 20 ml) connected to the existing Luer-Lock while keeping the three suction interrupts closed with a finger.

Clean the suction channel with a suitable cleaning brush (see TAB1-D-H) under running water, making sure that the openings at the ends are unobstructed. The tip of the cleaning brush must come out the end of the suction channel.

Flush the suction tube again five times with a disposable syringe (minimum volume 20 ml) connected to the existing Luer-Lock while keeping the three suction interrupts closed with a finger.

Machine cleaning and disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM (“Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.”, German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

- Ensure that the selected program is suitable for the instruments and includes sufficient flushing cycles,
 - that only sterile or low-germ (max. 1 germs/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (for instance purified water / highly purified water) is used for the final rinse,
 - the air used for drying is filtered, and
 - the disinfectior is inspected and maintained regularly.
- In selecting the cleaning agent system to be used, ensure that
- it is suitable for cleaning instruments made of metals and plastics,
 - a suitable disinfectant with proven effectiveness (for instance VAH (“Verbund für Angewandte Hygiene e.V.”, Applied Hygiene Association)/DGHM (“Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.”, German Association for Hygiene and Microbiology) approval, CE marking) is used and this is compatible with the chosen cleaning agent, and
 - the chemicals being used are compatible with the instruments (see “Material resistance”).

⚠ Complying with the concentrations, exposure times, flushing instructions and so on provided by the manufacturer of the cleaning agent or disinfectant is mandatory.

Process:

- Load the instruments into the disinfectior. In doing so, make sure the instruments do not touch each other. Using the Luer-Lock connection provided, connect the suction tube to the flushing connection of the disinfectior. Verify that openings are unobstructed prior to flushing. Use a cleaning wire / brush if needed (see TAB1-D-H).
- Start the program. Note the validation parameters for the cleaning / disinfection cycle listed below.
- After the program ends, disconnect the instruments from the connections and take them out of the disinfectior.
- Inspect the instruments as soon as possible after removal and, if required, package them after a subsequent drying cycle in a clean place (also see “Inspection”).

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.5 % (5 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 1 minute
Thermal disinfection	90±2 °C, 5 minutes

Proof of fundamental suitability of the instruments for effective machine cleaning and disinfection was provided by an independent, accredited test laboratory using the G 7836 CD disinfectior (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent deconex® 28 ALKA ONE-x (Borer Chemie AG, Zuchwil, Switzerland). Test report no. 121627-10 dated 25.04.2012.

Inspection:

Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

Maintenance:

None

Packaging:

Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packa-ging), or store in a suitable sterilization container meeting the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
- Sufficient protection of the instruments and/or sterilization packaging against mechanical damage
- Regular maintenance according to the information provided by the sterilization container manufacturer

Sterilization:

Products must be cleaned and disinfected prior to sterilization.

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	Min. 3 minutes
Drying time	30 minutes

Proof of fundamental suitability of the instruments for effective steam sterilization was provided by an inde-pendent, accredited test laboratory, test report no. 111739-10 dated 07.06.2011. Typical conditions in clinics and medical practices as well as the process described above were taken into account.

⚠ Do not sterilize in STERRAD®!

⚠ Do not sterilize in hot air!

⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

Material resistance:

Make sure that the substances listed below are not part of the cleaning agent or disinfectant:

- Organic acids, mineral acids, oxidizing acids (minimum allowable pH value 5.5)
- Strong lye (maximum allowable pH value 12)
- Organic solvents (such as acetone, ether, alcohol, gasoline)
- Oxidants (such as peroxide)
- Halogens (such as chlorine, iodine, bromine)
- Aromatic halogenated hydrocarbons

Do not use metal brushes, steel wool, scouring sponges or abrasive cleaning agents to clean instruments or sterilization trays and containers!

Do not expose instruments or sterilization trays and containers to temperatures in excess of 141 °C!


Storage / Transportation:

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging. For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medi-zintechnik.

Changes reserved. Current version available at www.sutter-med.de.

STERRAD® is a trademark of Johnson & Johnson, Inc.

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