REF 899000 | 2022-03-24 © Sutter Medizintechnik GmbH

Bipolare Pinzetten / Bipolar Forceps

Classic/Classic Micro

SELECTAL™

SuperGliss®/SuperGliss®ELP

SuperGliss®TEO/

SuperGliss®zhora

Irrigation/Suction

REF:

700150 – 700399 - incl. S, SV, SV1, SV2, F

700700 – 700799 - incl. S, SV


700800 – 700899 - incl. S, SV

702100 – 702299 - incl. S, SV

780130 – 786999 - incl. SG, SGS, SGSV,

SGSSV, SL, SLS, SGZ, SGSZ

TAB1		REF
A		992901012 g 290 mm, Ø 1.2 mm 992901018 g 290 mm, Ø 1.8 mm 993801018 g 380 mm, Ø 1.8 mm
B		701724, 701725 g n=2 701764, 701765 g n=5 701766, 701767 g n=10

DE	Bipolare Pinzetten	5
EN	Bipolar Forceps	9
FR	Forceps Bipolaires	13
ES	Pinzas Bipolares	17
IT	Pinze Bipolari	22
NL	Bipolaire Pincet	26
PT	Pinças Bipolar	30
SV	Bipolära Pincetter	35
DK	Bipolære Pincetter	39

FI	Bipolaariset Pinsetit	43
NO	Bipolære Pinsetter	47
PL	Bipolarne Pęsety	51
BG	Биполярни Пинсети	55
RO	Pense Bipolare	60
SK	Bipolárne Pinzety	65
SL	Bipolarne Pincete	69
CZ	Bipolární Pinzety	73
HU	Bipoláris Csipeszek	77

LT	Bipoliariniai Pincetai	81
LV	Bipolārās Pincetes	85
EL	Διπολική Λαβίδα	89
HR	Bipolarne pincete	93
RU	Биполярные Пинцеты	97
SR	Bipolarne pincete	102
MK	Биполарни пинцети	106
TR	Bipolar Pensle	111

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STERRAD® is a trademark of Johnson & Johnson, Inc.

Product / User / Disposal:

Electrosurgery accessories may only be used and disposed of by qualified medical personnel! Dispose of the instrument according to internal hospital guidelines for sharp, biologically contaminated items.

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory. For further information about electrical safety, we recommend IEC TR 61289 or DIN EN 60601-2-2 supplement 1.

⚠ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

Intended purpose:

Sutter SuperGliss®, SuperGliss®ELP, SuperGliss®TEO, SuperGliss®zhora, Classic, Classic Micro and Selectal™ bipolar forceps for coagulation of the selected tissue.

Sutter bipolar suction forceps for use in electrosurgery for the coagulation of tissue and suctioning fluids.

Sutter bipolar irrigation forceps for use in electrosurgery for coagulation and to supply fluid to select tissue.

Service life:

A minimum of 20 reconditioning cycles can be expected with proper use.

Before using:

⚠ Inspect the product for cleanliness, mechanical function, and intact insulation before each use.

We recommend checking the insulation with a suitable test device.

⚠ Only use sterilized products that are in flawless condition!

A certain discoloration of the non-stick instrument tip is normal and harmless.

To connect the forceps and cables, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

Electrosurgery cables:

The Sutter bipolar forceps are intended for use with bipolar silicone cables with US male connector or European flat connector from the manufacturer Sutter Medizintechnik GmbH.

During use:

⚠ Always work with the lowest output setting for the desired surgical effect.

⚠ Maximum allowable voltage 500 Vp.

⚠ Regularly wipe blood and tissue residue from the tips.

⚠ Tips of the forceps may cause injuries!

⚠ After application, the tips of the forceps may be so hot they cause burns!

⚠ Never lay down instruments on or in the immediate vicinity of the patient! Lay cables and store unused instruments so they are isolated from the patient.

⚠ Not for use in the presence of combustible or explosive substances!

Reconditioning:

General information:

Observe national guidelines and regulations!

Disconnect the instrument from the cable!

The overall reconditioning process encompasses pre-cleaning, cleaning / disinfection, and sterilization.

⚠ Machine cleaning / disinfection is preferred due to effectiveness and reproducibility!

⚠ Irrigation and suction forceps must always be reconditioned by machine!

⚠ Do not immerse in hydrogen peroxide (H₂O₂)!

⚠ Do not bend open the forceps! **(FIG1)**

Pre-cleaning:

- Do not allow blood and tissue residues to dry up; rinse thoroughly with cold water after max. 1 hour! Use a soft brush if needed (no wire brush or similar).

- On irrigation and suction forceps, verify that openings are unobstructed prior to flushing. Use a cleaning wire / brush if needed (**TAB1:A**). Flush the lumina of the instruments five times with a disposable syringe (minimum volume 50 ml) and direct connection to the existing Luer-Lock.
- Moveable parts must be moved back and forth several times during pre-cleaning.

Manual cleaning and disinfection:

Cleaning step	Description
Pre-cleaning	Flush for 5 minutes under cold water, operating moving parts. Clean the instrument with a soft brush (such as MED100.33 from Medisafe GmbH) until no more residues are visible.
Ultrasound and disinfection	Ultrasound bath 35 kHz at room temperature, 10 minutes, cleaning / disinfection solution 2 % Bomix® plus (Bode Chemie).
Secondary cleaning	Rinse hard to reach areas for 20 seconds with a spray nozzle as needed, then rinse the entire instrument for 30 seconds with demineralized water.

Machine cleaning and disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

- Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported.
- The optionally available storage trays (**TAB1:B**) ensures secure support.
For improved cleaning results, the lumina of the irrigation and suction forceps can be connected to the flushing connection of the CDD using the existing Luer-Lock connections.

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.5 % (5 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 1 minute
Thermal disinfection	90±2 °C, 5 minutes

- Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the instruments and are permitted according to the A₀ concept, comparable A₀ value > 3000. If a different cleaning agent is used, only use products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie), for instance regarding the pH value and compatibility with plastics. Contact the responsible supplier or hygiene officer in case of doubt.

Inspection:

Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

Maintenance:

None

Packaging:

Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packaging), or wrap the instrument or tray with the cleaned and disinfected instruments in a cotton cloth and store them together in suitable sterilization containers meeting the following requirements:

- EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
- Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

Sterilization:

Products must be cleaned and disinfected prior to sterilization.

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	min. 3 minutes
Drying time	min. 30 minutes

⚠ Do not sterilize in hot air!

⚠ Do not sterilize in STERRAD®!

⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

Storage / Transportation:

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging.

For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Please note:

The instructions listed above have been validated by the manufacturer as suitable for preparing a medical device for reuse. It is the responsibility of the reconditioner to ensure that actual reconditioning with the equipment, materials, and personnel used in the reconditioning facility achieves the desired result.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The product may only be repaired by the manufacturer or an agent expressly authorized by the manufacturer. Otherwise the warranty and any possible additional liability claims against the manufacturer are voided.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik GmbH.

Changes reserved. Current version available at www.sutter-med.de.

STERRAD® is a trademark of Johnson & Johnson, Inc.

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