REF 899002 | 2021/07 © Sutter Medizintechnik GmbH

# Bipolare Elektroden / Bipolar Electrodes

REF:  
700462 – 700499

<b>DE</b> Bipolare Elektroden	5	<b>FI</b> Bipolaariset Elektrodit	39
<b>EN</b> Bipolar Electrodes	9	<b>NO</b> Bipolære Elektroder	43
<b>FR</b> Électrodes Bipolaires	12	<b>PL</b> Elektrody Bipolarne	46
<b>ES</b> Electrodos Bipolares	16	<b>BG</b> Биполярни Електроди	50
<b>IT</b> Elettrodi Bipolari	21	<b>RO</b> Electrozi Bipolari	54
<b>NL</b> Bipolaire Elektroden	24	<b>SK</b> Bipolárne Elektródy	58
<b>PT</b> Elétrodos Bipolares	28	<b>SL</b> Bipolarne Elektrode	62
<b>SV</b> Bipolär Blektrod	33	<b>CZ</b> Bipolární Elektrody	65
<b>DK</b> Bipolære Blektroder	36	<b>HU</b> Bipoláris Elektródák	69

<b>LT</b>	Bipoliariniai Elektrodai	73
<b>LV</b>	Bipolārie Elektrodi	76
<b>EL</b>	Διπολικά Ηλεκτρόδια	80
<b>HR</b>	Bipolarne elektrode	84
<b>RU</b>	Биполярные Электроды	87
<b>SR</b>	Bipolarne elektrode	91
<b>MK</b>	Биполарни електроди	95
<b>TR</b>	Bipolar Elektrotlar	99

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STERRAD® is a trademark of Johnson & Johnson, Inc.

## Product / User / Disposal:

Electrosurgery accessories may only be used and disposed of by qualified medical personnel! Dispose of the instrument according to internal hospital guidelines for sharp, biologically contaminated items.

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.

⚠ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

## Intended purpose:

Bipolar electrodes for the coagulation / volume reduction / cutting of soft tissue during surgical procedures.

## Service life:

A minimum of 20 reconditioning cycles can be expected with proper use.

## Before using:

⚠ Inspect the product for cleanliness, intact insulation, and integrity before each use. Note the sensitive instrument tips in particular.

⚠ Only use sterilized products that are in flawless condition!

All products must be used with a suitable connection cable from Sutter!

⚠ To connect electrodes or accessories, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

## During use:

Operating guidelines for these product groups are available as non-binding instructions. They can be requested from the manufacturer.

⚠ Always work with the lowest output setting for the desired surgical effect.

⚠ Electrode tips may cause injuries!

⚠ After application, electrode tips may be so hot they cause burns!

⚠ Never lay down the electrode on or in the immediate vicinity of the patient! Lay cables and store unused instruments so they are isolated from the patient.

⚠ Not for use in the presence of combustible or explosive substances!

⚠ Unintentional activation or movement of the activated electrode outside the field of vision can lead to unwanted burns!

⚠ Maximum allowable voltage 500 Vp.

## Reconditioning:

Observe national guidelines and regulations!

Disconnect electrodes from accessories / connection cables!

The overall reconditioning process encompasses pre-cleaning, cleaning / disinfection, and sterilization.

⚠ Machine cleaning / disinfection is preferred due to effectiveness and reproducibility!

⚠ Do not immerse in hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>)!

#### Pre-cleaning:

Do not allow blood and tissue residues to dry up; rinse thoroughly with cold water after max. 1 hour! Use a soft brush if needed (no wire brush or similar).

#### Manual cleaning and disinfection:

Cleaning step	Description
Pre-cleaning	Rinse for 5 minutes under cold water and clean with a soft brush (such as MED100.33 from Medisafe GmbH) until no more residues are visible.
Ultrasound and disinfection	Ultrasound bath 35 kHz at room temperature, 10 minutes, cleaning/ disinfection solution 2 % Bomix® plus (Bode Chemie).
Secondary cleaning	Rinse hard to reach areas for 20 seconds with a spray nozzle as needed, then rinse the entire instrument for 30 seconds with demineralized water.

#### Machine cleaning and disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have

been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

- Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported.

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.3 % (3 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 2 minutes
Thermal disinfection	90±2 °C, 5 minutes

- Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the instruments and are permitted according to the A<sub>0</sub> concept, comparable A<sub>0</sub>-value>3000. If a different cleaning agent is used, only use products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie), for instance regarding the pH value and compatibility with plastics. Please contact your supplier or hygiene officer in case of questions.

**Inspection:**

Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

**Maintenance:**

None

**Packaging:**

Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packaging – note sharp instrument tips!), or store the instrument in a suitable sterilization container meeting the following requirements:

- EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
- Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

**Sterilization:**

Products must be cleaned and disinfected prior to sterilization.

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	min. 3 minutes
Drying time	min. 30 minutes

⚠ Do not sterilize in hot air!

⚠ Do not sterilize in STERRAD®!

⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

**Storage / Transportation:**

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging. For returns, products must be cleaned, disinfected, and packed in sterile packaging.

**Please note:**

The instructions listed above have been validated by the manufacturer as suitable for preparing a medical device for reuse. It is the responsibility of the reconditioner to ensure that actual reconditioning with the equipment, materials, and personnel used in the reconditioning facility achieves the desired result.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The product may only be repaired by the manufacturer or an agent expressly authorized by the manufacturer. Otherwise the warranty and any possible additional liability claims against the manufacturer are voided.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik GmbH.

Changes reserved. Current version available at [www.sutter-med.de](http://www.sutter-med.de).

STERRAD® is a trademark of Johnson & Johnson, Inc.

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