Monopolare Elektroden / Monopolar Electrodes


Standard
ARROWtip™
Arthro
Micro

REF:

360320 – 360379
360342 SV1
360440 – 360463
360500 – 360529
360531

360535
360540 – 360549
360800 – 360818
360830 – 36085

TAB1				
	REF	Max. Vp	Ø [mm]	Ø [mm]
	360320 – 360379 360342 SV1	3000	2,4 (3/32")	2,4 (3/32")
	360500 – 360521	1000	1,6	1,6
	360531 / 360535	900	1,6	1,6 Sutter-REF 360215
	360540 – 360549	2000	1,6	1,6
	360800 – 360818	2000	2,4 (3/32")	2,4 (3/32")
	360830 – 360859	500	2,4 (3/32")	2,4 (3/32")
	360885 – 360893	2000	2,4 (3/32")	2,4 (3/32")

DE	Monopolare Elektroden	5
EN	Monopolar Electrodes	9
FR	Électrodes Monopolaires	13
ES	Electrodos Monopolares	17
IT	Elettrodi Monopolari	20
NL	Monopolaire Elektroden	24
PT	Elétrodos Monopolares	28
SV	Monopolära Elektroder	32
DK	Monopolære Elektroder	36

FI	Monopolaariset Elektrodit	39
NO	Monopolære Elektroder	43
PL	Elektrody Monopolarne	46
BG	Монополярни Електроди	50
RO	Electrozi Monopolari	54
SK	Monopolárne Elektródy	58
SL	Monopolarne Elektrode	62
CZ	Monopolární Elektrody	65
HU	Monopoláris Elektródák	69

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STERRAD® is a trademark of Johnson & Johnson, Inc.

Product / User / Disposal:

Electrosurgery accessories may only be used and disposed of by qualified medical personnel! Dispose of the instrument according to internal hospital guidelines for sharp, biologically contaminated items.

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory. For further information about electrical safety, we recommend IEC TR 61289 or DIN EN 60601-2-2 supplement 1.

⚠ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

Intended purpose:

Monopolar electrode for cutting and coagulating soft tissue during surgical procedures.

Service life:

The service life may be non-specifically shortened depending on the type of application. Replace the electrode after no more than 20 applications.

ARROWtip™ models feature dual-layer insulation. If the outer, black insulation is damaged and the inner, yellow insulation is visible, dispose of the electrode for safety reasons!

Before using:

⚠ Inspect the product for cleanliness, intact insulation, and damage before each use.

For long electrodes, we recommend checking the insulation with a suitable test device.

⚠ Only use sterilized products that are in flawless condition!

Insert the electrode connection fully and properly into the electrosurgery handle intended for the purpose.

To connect the electrode, the electrosurgery device must be turned off or in standby mode.

⚠ Failure to comply may lead to burns and electric shocks!

Verify that the neutral electrode is applied correctly!

During use:

⚠ Always work with the lowest output setting for the desired surgical effect.

⚠ Note the maximum voltages and connection information (**TAB1**).

⚠ Electrode tips may cause injuries!

⚠ After application, electrode tips may be so hot they cause burns!

⚠ Never lay down the electrode and/or handle on or in the immediate vicinity of the patient! Lay cables and store unused instruments so they are isolated from the patient.

⚠ Not for use in the presence of combustible or explosive substances!

When flushing instruments are used simultaneously with the electrode, use non-conductive flushing fluid if possible!

Reconditioning:

Observe national guidelines and regulations!

Disconnect the electrode / adapter / handle / cable from each other.

The overall reconditioning process encompasses pre-cleaning, cleaning / disinfection, and sterilization.

⚠ Machine cleaning / disinfection is preferred due to effectiveness and reproducibility!

⚠ Do not immerse in hydrogen peroxide (H₂O₂)!

Pre-cleaning:

Do not allow blood and tissue residues to dry up; rinse thoroughly with cold water after max. 1 hour! Use a soft brush if needed (no wire brush or similar).

Manual cleaning and disinfection:

Cleaning step	Description
Pre-cleaning	Rinse for 5 minutes under cold water and clean with a soft brush (such as MED100.33 from Medisafe GmbH) until no more residues are visible.
Ultrasound and disinfection	Ultrasound bath 35 kHz at room temperature, 10 minutes, cleaning/disinfection solution 2 % Bomix® plus (Bode Chemie).

Cleaning step	Description
Secondary cleaning	Rinse hard to reach areas for 20 seconds with a spray nozzle as needed, then rinse the entire instrument for 30 seconds with demineralized water.

Machine cleaning and disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

- Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported.

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.5 % (5 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 1 minute
Thermal disinfection	90±2 °C, 5 minutes

- Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the instruments and are permitted according to the A₀ concept, comparable A₀ value > 3000). If a different cleaning agent is used, only use products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie), for instance regarding the pH value and compatibility with plastics. Contact the responsible supplier or hygiene officer in case of doubt.

Inspection:

Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

Maintenance:

None

Packaging:

Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packaging – note sharp instrument tips!), or store the instrument in a suitable sterilization container meeting the following requirements:

- EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
- Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

Sterilization:

Products must be cleaned and disinfected prior to sterilization.

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	Min. 3 minutes
Drying time	Min. 15 minutes

⚠ Do not sterilize in hot air!

⚠ Do not sterilize in STERRAD®!

⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

Storage / Transportation:

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging.
For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Please note:

The instructions listed above have been validated by the manufacturer as suitable for preparing a medical device for reuse. It is the responsibility of the reconditioner to ensure that actual reconditioning with the equipment, materials, and personnel used in the reconditioning facility achieves the desired result.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The product may only be repaired by the manufacturer or an agent expressly authorized by the manufacturer. Otherwise the warranty and any possible additional liability claims against the manufacturer are voided.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizin-technik GmbH.

Changes reserved. Current version available at www.sutter-med.de.

STERRAD® is a trademark of Johnson & Johnson, Inc.

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Sutter Medizintechnik GmbH