REF 899010 | 2021/09 © Sutter Medizintechnik GmbH

# Monopolare Saugrohre / Monopolar Suction Tubes

REF:  
715000 – 715019

TAB 1			REF
A		Sutter (BM-780 II), Martin, Berchtold, Aesculap, Erbe (T-Serie)	360185
B		Erbe ACC / IC / ViO, Karl Storz	360186
C		Sutter CURIS®, Valleylab, Conmed, Bovie, Bowa	360187
D		290 mm, Ø 3.2 mm -->  (REF 715010)	992901032
E		290 mm, Ø 2.2 mm -->  (REF 715015)	992901022
F		380 mm, Ø 2.2 mm -->  (REF 715017)	993801022
G		290 mm, Ø 3.2 mm -->  (REF 715018)	992901032
H		290 mm, Ø 3.2 mm -->  (REF 715019)	992901032

<b>DE</b> Monopolare Saugrohre	5	<b>FI</b> Monopolaariset Imuputket	51
<b>EN</b> Monopolar Suction Tubes	10	<b>NO</b> Monopolære Sugerør	56
<b>FR</b> Tuyau D'aspiration Monopolaire	15	<b>PL</b> Monopolarne Rury Ssawne	61
<b>ES</b> Tubos Monopolares de Aspiración	20	<b>BG</b> Монополярни Смукателни Тръби	66
<b>IT</b> Tubi di Aspirazione Monopolari	25	<b>RO</b> Tuburi Aspirație Monopolare	72
<b>NL</b> Monopolaire Zuigbuizen	30	<b>SK</b> Monopolárne Sacie Rúry	77
<b>PT</b> Tubos de Aspiração Monopolares	36	<b>SL</b> Monopolarne Sesalne Cevi	82
<b>SV</b> Monopolära Sugrör	41	<b>CZ</b> Monopolární Sací Trubičky	87
<b>DK</b> Monopolære Sugerør	46	<b>HU</b> Monopoláris Szívócsövek	92

<b>LT</b>	Monopoliariniai Siurbimo Vamzdeliai	97
<b>LV</b>	Monopolārās Sūcējcaurules	102
<b>EL</b>	Μονοπολικοί Σωλήνες Αναρρόφησης	107
<b>HR</b>	Monopolarne usisne cijevi	112
<b>SR</b>	Monopolarne usisne cevi	117
<b>MK</b>	Μονοποларни цевки за вшмукување	122
<b>TR</b>	Monopolar Emme Boruları	128
<b>ZH</b>	单极吸引管	133
<b>ID</b>	Tabung Isap Monopolar	137

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STERRAD® is a trademark of Johnson & Johnson, Inc.

**Product / User / Disposal:**

Electrosurgery accessories may only be used and disposed of by qualified medical personnel! Dispose of the instrument according to internal hospital guidelines for sharp, biologically contaminated items.

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory. For further information about electrical safety, we recommend IEC TR 61289 or DIN EN 60601-2-2 supplement 1.

⚠ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

**Intended purpose:**

Monopolar suction tube for suctioning smoke and fluids and for coagulating soft tissue during surgical procedures.

**Service life:**

A minimum of 20 reconditioning cycles can be expected with proper use.

**Before using:**

⚠ Inspect the product for cleanliness, intact insulation, and integrity before each use.

We recommend checking the insulation with a suitable test device.

⚠ Only use sterilized products that are in flawless condition!

A certain discoloration of the tips on non-stick instruments is normal and harmless.

Do not use in case of discernible damage to the electrical insulation, visible corrosion, or if the product exhibits nicks, dents, chipping, or cracks. Properly dispose of the product in such cases.

⚠ Only use the instrument with a suitable connection cable from Sutter! (see **TAB1:A-C**).

To connect the product, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

⚠ Never use monopolar instruments without a neutral electrode! Verify that the neutral electrode is applied correctly.

⚠ Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory. Carefully follow the instructions provided by the manufacturer of the electrosurgery device being used, and read all warnings and safety notices.

**During use:**

⚠ Always work with the lowest output setting for the desired surgical effect.

Operation only in COAG mode with max. 40 watts. Brief, pulse-like activation if possible. Avoid continuous activation > 15 seconds.

Do not use in spray / fulguration mode!

Avoid spark formation!

- ⚠ After application, the tips may be so hot they cause burns!
- ⚠ Never lay down instruments on or in the immediate vicinity of the patient!
- ⚠ Position cables so the possibility of contact with the patient and other electrical lines is excluded.
- ⚠ When instruments are temporarily unused, store them so they are isolated from the patient.
- ⚠ Not for use with or in the vicinity of combustible or explosive substances!
- ⚠ When flushing instruments are used simultaneously, use non-conductive flushing fluid if possible!
- ⚠ Maximum allowable voltage 3000 Vp.

#### **Reconditioning:**

##### **General information:**

Observe national guidelines and regulations!

Disconnect the instrument from the cable!

Do not allow blood and tissue residues to dry up!

Remove blood and tissue residues with a soft cloth or brush!

Do not use sharp tools / scouring agents!

⚠ Do not immerse in hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>)!

##### **Manual cleaning and disinfection:**

⚠ The instrument must always be reconditioned by machine – no manual cleaning!

According to the recommendations of the DGSV ("Deutsche Gesellschaft für Sterilgutversorgung", German Association for Sterile Goods Supply), monopolar suction tubes are assigned to risk group B\*. Machine cleaning is required for such products.

*\*This classification was assigned according to the DGSV 2013 flow chart for the classification of medical devices based on the recommendations of the KRINKO ("Krankenhaushygiene und Infektionsprävention Kommission", Hospital Hygiene and Infection Prevention Commission)/BfArM ("Bundesinstitut für Arzneimittel und Medizinprodukte", Federal Institute for Drugs and Medical Devices), Bundesgesundheitsblatt (Federal Health Bulletin) 2012; 55:1244-1310*

##### **Manual pre-cleaning:**

To make subsequent cleaning easier, do not allow tissue residue and blood to dry on during use.

Clean the instruments to remove major contamination directly after use (within 1 hour at most). Use running water or a disinfectant solution. The disinfectant must be aldehyde-free (to prevent fixing blood particles) with proven effectiveness (for instance VAH ("Verbund für Angewandte Hygiene e.V", Applied Hygiene Association)/DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) approval, CE marking) and be suitable for disinfecting the instruments (see "Material resistance").

Only use a soft brush or soft, clean cloth to manually remove contaminants – no metal brushes, steel wool, scouring sponges, or abrasive cleaning agents.

Flush the suction tube five times with a disposable syringe (minimum volume 20 ml) connected to the existing Luer-Lock while keeping the three suction interrupts closed with a finger.

Clean the suction channel with a suitable cleaning brush (see **TAB1:D-H**) under running water, making sure that the openings at the ends are unobstructed. The tip of the cleaning brush must come out the end of the suction channel.

Flush the suction tube again five times with a disposable syringe (minimum volume 20 ml) connected to the existing Luer-Lock while keeping the three suction interrupts closed with a finger.

#### **Machine cleaning and disinfection:**

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

- Ensure that the selected program is suitable for the instruments and includes sufficient flushing cycles,
- that only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (for instance purified water / highly purified water) is used for the final rinse,
- the air used for drying is filtered, and
- the disinfectant is inspected and maintained regularly.

In selecting the cleaning agent system to be used, ensure that

- it is suitable for cleaning instruments made of metals and plastics,

- a suitable disinfectant with proven effectiveness (for instance VAH ("Verbund für Angewandte Hygiene e.V", Applied Hygiene Association)/DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) approval, CE marking) is used and this is compatible with the chosen cleaning agent, and
- the chemicals being used are compatible with the instruments (see "Material resistance").

⚠ Complying with the concentrations, exposure times, flushing instructions and so on provided by the manufacturer of the cleaning agent or disinfectant is mandatory.

Process:

1. Load the instruments into the disinfectant. In doing so, make sure the instruments do not touch each other. Using the Luer-Lock connection provided, connect the suction tube to the flushing connection of the disinfectant. Verify that openings are unobstructed prior to flushing. Use a cleaning wire / brush if needed (see **TAB1:D-H**).
2. Start the program. Note the validation parameters for the cleaning / disinfection cycle listed below.
3. After the program ends, disconnect the instruments from the connections and take them out of the disinfectant.
4. Inspect the instruments as soon as possible after removal and, if required, package them after a subsequent drying cycle in a clean place (also see "Inspection").

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.5 % (5 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 1 minute
Thermal disinfection	90±2 °C, 5 minutes

**Inspection:**

Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

**Maintenance:**

None

**Packaging:**

Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packaging), or store in a suitable sterilization container meeting the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
- Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

- Regular maintenance according to the information provided by the sterilization container manufacturer

**Sterilization:**

Products must be cleaned and disinfected prior to sterilization.

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	Min. 3 minutes
Drying time	Min. 30 minutes

⚠ Do not sterilize in STERRAD®!

⚠ Do not sterilize in hot air!

⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

**Material resistance:**

Make sure that the substances listed below are not part of the cleaning agent or disinfectant:

- Organic acids, mineral acids, oxidizing acids (minimum allowable pH value 5.5)
- Strong lye (maximum allowable pH value 12)
- Organic solvents (such as acetone, ether, alcohol, gasoline)
- Oxidants (such as peroxide)
- Halogens (such as chlorine, iodine, bromine)
- Aromatic halogenated hydrocarbons

Do not use metal brushes, steel wool, scouring sponges or abrasive cleaning agents to clean instruments or sterilization trays and containers!

Do not expose instruments or sterilization trays and containers to temperatures in excess of 141 °C!

**Storage / Transportation:**

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging.

For returns, products must be cleaned, disinfected, and packed in sterile packaging.

**Please note:**

The instructions listed above have been validated by the manufacturer as suitable for preparing a medical device for reuse. It is the responsibility of the reconditioner to ensure that actual reconditioning with the equipment, materials, and personnel used in the reconditioning facility achieves the desired result.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The product may only be repaired by the manufacturer or an agent expressly authorized by the manufacturer. Otherwise the warranty and any possible additional liability claims against the manufacturer are voided.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik GmbH.

Changes reserved. Current version available at [www.sutter-med.de](http://www.sutter-med.de).

STERRAD® is a trademark of Johnson & Johnson, Inc.

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