BM-780 II

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The BM-780 II radiofrequency generator is ideally suited for small surgical interventions in ENT clinics and doctors' offices, for plastic surgery, and dermatology.

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AutoRF™ is a smart impedance control function that will tailor the RF power output to the tissue condition. Whether it is cutting through different types of tissue or altering tissue conditions during coagulation, the AutoRF™ feature will deliver adapted power output as required by the different tissue impedances.

Thanks to AutoRF™ the BM-780 II radiofrequency generator offers more safety and convenience. Cutting and coagulation results are reproducible.

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

In bipolar coagulation the power output is reduced according to the degree of coagulation and the desiccation of tissue to a minimum of only about ¼ of the initial setting.

During the monopolar cutting process the unit adjusts the output to the type of tissue (skin, muscles, fat, etc.) as needed.

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„R“ signifies electric tissue resistance and „P“ the actual power output. Sections 1 to 5 show the different kinds of tissues and cutting speeds to which the unit adjusts its power output automatically.

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
What is RaVoR™ (Radiofrequency Volume Reduction)?

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Radiofrequency energy causes deliberate local lesions through the short and selective heating of certain tissue areas. The treated tissue is decomposed and transformed into fibrous scar tissue. This process leads to a shrinkage and stiffening of the treated area.

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RaVoR™ Bipolar Electrodes

The specifically designed bipolar electrodes for RaVoR™ ensure a lesion suited to the patient's anatomy. A special, ultra-thin insulation protects surface of the mucosa from thermal damage.

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RaVoR™

bipolar electrode for the inferior turbinate, single-use
working length: 103 mm


37 38 51

bipolarer Adapter Sutter BM-780 II
für RaVoR™ bipolare Elektroden, single-use
bipolar adapter Sutter BM-780 II
for RaVoR™ bipolar electrodes, single-use

“RaVoR™ is a modern surgical technique showing good and long-lasting treatment results when used to reduce the volume of hypertrophic turbinates. At the same time it preserves the mucosa and its function.”

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
Inferior turbinate – Preoperative condition and six months postoperatively with significantly enlarged nasal passage.

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Schematic view of the puncture sites for the application of radiofrequency energy of the hypertrophic turbinates.

RaVoR™ am weichen Gaumen | RaVoR™ of the soft palate

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RaVoR™

bipolar electrode for the soft palate, single-use
working length: 110 mm


37 38 51

bipolarer Adapter Sutter BM-780 II
für RaVoR™ bipolare Elektroden, single-use
bipolar adapter Sutter BM-780 II
for RaVoR™ bipolar electrodes, single-use

“The radiofrequency assisted soft palate procedure is a minimally invasive, safe and quick procedure. It is well tolerated by patients. We have not observed any bleeding that needed special attention.”

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
Puncture sites for the application of radiofrequency energy in the soft palate.

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Dissection of surplus uvula tissue and incision lines for the triangular excision of the posterior palatal pillars. (With ARROWtip™ monopolar microdissection electrode, REF: 36 44 42).

+

37 38 51
bipolarer Adapter Sutter BM-780 II für RaVoR™ bipolare Elektroden, single-use
bipolar adapter Sutter BM-780 II for RaVoR™ bipolar electrodes, single-use


70 44 99

RaVoR™
bipolare Elektrode für den Zungengrund, single-use
Arbeitslänge: 110 mm
RaVoR™
bipolar electrode for the tongue base, single-use
working length: 110 mm

“In our clinical practice we successfully make use of the tongue base radiofrequency procedure. Our experience shows that when using this minimally invasive method together with other surgical techniques, the outcome of sleep-related breathing disorder surgery can be improved. The treatment is useful and should be considered in the treatment of patients with tongue base collapse.”

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
Puncture sites for the treatment of the tongue base.

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The low profile of the instrument and its strong shaft enable the surgeon to insert the electrode at the back of the tongue.

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

RaVoR™

bipolar electrode for the posterior pillars
working length: 119 mm


Hinteres Gaumensegel unmittelbar nach Radiofrequenz-Behandlung.
Pillar immediately after radiofrequency treatment.

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Histology of the pillar: Coagulation and necrosis (a), thrombosis (b), atrophic muscle (c).

Gestrafftes Gaumensegel nach einer Woche postoperativ.
Retracted pillar one week post-operatively.

"For the posterior pillars alone two sessions may be required. It is noteworthy that there have not been any perioperative or post-operative complications and hardly any bleedings. This method is ideal for the office-based setting."

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RaVoR™

bipolar electrode for the tongue base tonsils
working length: 142 mm

Ideale Produktkombination für die Blutstillung der Adenoidektomie | Ideal product combination for coagulation of the adenoids

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***Disclaimer:** The information presented herein has been carefully researched and compiled with the help of specialist physicians. They are not meant to serve as a detailed treatment guide. They do not replace the user instructions for the medical devices used. Sutter accepts no liability for the treatment results beyond the mandatory legal regulations.

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Product availability is subject to regulatory approval in individual markets. Products may therefore not be available in all markets.
The listed lengths and sizes serve as a guideline and may be rounded up or down. The actual lengths may vary slightly.

PRECISION ELECTROSURGERY
Made in Germany

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ALFRED-WALZ-STR. 22 · 79312 EMMENDINGEN/GERMANY
TEL. +49(0)7641-96256-0 · FAX +49(0)7641-96256-30
WWW.SUTTER-MED.COM · INFO@SUTTER-MED.DE