**REF:**

37 01 00 R, 37 38 50, 37 50 30

36 01 48

37 01 30 A, G, K, L, P, R, S, T, W

37 01 32 L

37 01 35 A, G, K, L, P, R, S, T, W, SG

37 01 38 A, G, K, L, P, R, S, T, W

37 01 40 L, S, SG

37 01 45 G, L

37 01 50 A, G, K, L, P, R, S, T, VS, W, SG

37 01 53 A, G, K, L, P, R, S, T, W, SG

37 01 54 A, G, L, P, R, S, SG

36 85 87, 36 86 87

36 01 70...36 01 73

36 01 73-18

36 01 73-30

36 01 85...36 01 88

36 01 90...36 01 93

36 01 97

36 02 36...36 02 38

36 02 14, 36 02 15, 36 02 18

36 07 01...36 07 05

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STERRAD® is a trademark of Johnson & Johnson, Inc.

**English****Product / User / Disposal:**

Electrosurgery accessories may only be used and disposed of by qualified medical personnel! Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.

⚠ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

Intended use:

For the connection of electrosurgery accessories (such as electrodes, forceps) to a compatible electrosurgery device.

Service life:

A minimum of 20 reconditioning cycles can be expected with proper use.

Before using:

- ⚠ Inspect the product for cleanliness, intact insulation and integrity before each use.
- ⚠ Only use sterilized products that are in flawless condition!
- ⚠ Only for use with compatible products and plugs. Do not connect to different or unknown power sources!
- ⚠ Use the correct connection for the product on the electrosurgery device or instrument: Monopolar – bipolar – neutral electrode.
- ⚠ To connect the product, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

During use:

- ⚠ To plug in and unplug cables, grasp them only on the plug part. Do not kink cables. Cables may not be wound up or rolled up during use.
- Cables must not be improperly attached.
- ⚠ Not for use in the presence of combustible or explosive substances!
- ⚠ Note maximum allowable voltages: Monopolar cables: 5000 Vp Bipolar cables: 1000 Vp Monopolar handles: 4000 Vp

Reconditioning:**General information:**

Observe national guidelines and regulations!
Disconnect cables / handles from the instrument / electrosurgery device.
Do not use sharp tools / scouring agents for cleaning!
Clean the cables / handles under running cold water to remove major contamination directly after use (within 1 hour at most).

Machine cleaning / Disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883). During machine cleaning / disinfection, store the products securely and protected against mechanical damage. Do not clean together with sharp-edged or pointy objects.

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.5 % (5 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 1 minute
Thermal disinfection	90±2 °C, 5 minutes

Inspection:

Prior to sterilization, visually inspect the cable / handle and check for intact insulation, cleanliness, and integrity.

Maintenance:

None

Packaging:

- Package cleaned and disinfected cables / handles in disposable sterilization packaging or sterilization containers (single or double packaging) meeting the following requirements:
 - EN ISO/ANSI AAMI ISO 11607
 - Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
 - Sufficient protection of the cable / handles and/or sterilization packaging against mechanical damage

Sterilization:

Products must be cleaned and disinfected prior to sterilization.

Use **only** the sterilization processes described in the following:

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated pre-vacuum process
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665. Max. sterilization time 20 minutes)
Sterilization time (holding time at sterilization temperature)	4 minutes

Proof of fundamental suitability of the instruments for effective steam sterilization was provided by an independent, accredited test laboratory, test report no. 121741 3776 dated 13.12.2002. Typical conditions in clinics and medical practices as well as the process described above were taken into account.

Always use the following sterilization process for the 36 01 85/86/87 cables:

Program steps	Parameters
Process	Steam sterilization
Sterilization temperature	134 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	20 minutes

Proof of fundamental suitability of the instruments for effective steam sterilization was provided by an independent, accredited test laboratory, test report no. 99d005 dated 26.04.1999. Typical conditions in clinics and medical practices as well as the process described above were taken into account.

Ensure adequate drying!

- ⚠ Do not sterilize in hot air!
- ⚠ Do not sterilize in STERRAD®!
- ⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the products and do not use them again.

Storage / Transportation:

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging.

Roll cables up loosely; do not wind tightly, kink, or fold.

For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Please note:

The 37 01 35L/K/S/A/G/R/T/W/SG cables do not correspond to section 201.12.2f of IEC 60601-2-2:2009.

- ⚠ Only connect these cables to the bipolar output of electrosurgery devices. The risk is reasonable given proper use by qualified personnel.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik.

Changes reserved, current version available under www.sutter-med.de.

STERRAD® is a trademark of Johnson & Johnson, Inc.

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