**Classic/Classic Micro SELECTAL™ SuperGliss®/SuperGliss®ELP SuperGliss®TEO/SuperGliss®zhora**  
**Irrigation/Suction**

**REF:**

- 70 01 50 – 70 03 99 – incl. S, SV, SV1, SV2, F
- 70 07 00 – 70 07 99 – incl. S, SV
- 70 08 00 – 70 08 99 – incl. S, SV
- 70 21 00 – 70 22 99 – incl. S, SV
- 78 01 30 – 78 69 99 – incl. SG, SGS, SGSV, SGSSV, SL, SLS, SGZ, SGSZ

TAB 1		REF
A		992901012, 992901018, 993801018
B		701724, 701725, 701764, 701765, 701766, 701767
C		99 00 80
D		80 00 02


FIG. 1

**Deutsch**

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**English**

**Product / User / Disposal:**  
 Electroscopy accessories may only be used and disposed of by qualified medical personnel! Notwithstanding these instructions, reading the instructions for use provided with the electroscopy device and other accessories being used is mandatory.

△ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

**Intended use:**  
 Sutter SuperGliss®, SuperGliss®ELP, SuperGliss®TEO, SuperGliss®zhora, Classic, Classic Micro and Selectal™ bipolar forceps for coagulation of the selected tissue.  
 Sutter bipolar suction forceps for use in electroscopy for the coagulation of tissue and suctioning fluids.  
 Sutter bipolar suction forceps for use in electroscopy for coagulation and to supply fluid to select tissue.

**Service life:**  
 A minimum of 20 reconditioning cycles can be expected with proper use.

**Before using:**  
 △ Inspect the product for cleanliness, mechanical function, and intact insulation before each use.  
 We recommend checking the insulation with a suitable test device.  
 △ Only use sterilized products that are in flawless condition!  
 A certain discoloration of the non-stick instrument tip is normal and harmless.  
 To connect the forceps and cables, the electroscopy device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!  
 For further information about electrical safety, we recommend DIN EN 60601-2-2 supplement 1.

**Electrosurgery cables:**  
 The Sutter bipolar forceps are intended for use with bipolar silicone cables with US male connector or European flat connector from the manufacturer Sutter Medizintechnik GmbH.

**During use:**  
 Always work with the lowest output setting for the desired surgical effect.  
 △ Maximum allowable voltage 500 Vp.  
 △ Regularly wipe blood and tissue residue from the tips.  
 △ Tips of the forceps may cause injuries!  
 △ After application, the tips of the forceps may be so hot they cause burns!  
 △ Never lay down instruments on or in the immediate vicinity of the patient! Lay cables and store unused instruments so they are isolated from the patient.  
 △ Not for use in the presence of combustible or explosive substances!

**Reconditioning:**  
**General information:**  
 Observe national guidelines and regulations!  
 △ Verify the propreté, le fonctionnement mécanique et l'intégrité parfaites de l'appareil avant chaque utilisation.  
 The overall reconditioning process encompasses pre-cleaning, cleaning / disinfection, and sterilization.  
 △ Machine cleaning / disinfection is preferred due to effectiveness and reproductibility!  
 △ Flushing and suction forceps must always be reconditioned by machine!  
 △ Do not immerse in hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>)!  
 △ Do not bend open the forceps! (FIG1)

**Pre-cleaning:**  
 • Do not allow blood and tissue residues to dry up; rinse thoroughly with cold water after max. 1 hour! Use a soft brush if needed (no wire brush or similar).  
 • On flushing and suction forceps, verify that openings are unobstructed prior to flushing. Use a cleaning wire / brush if needed (TAB1-A). Flush the lumina of the instruments five times with a disposable syringe (minimum volume 50 ml) and direct connection to the existing Luer-Lock.  
 • Moveable parts must be moved back and forth several times during pre-cleaning.

Cleaning step	Description
Pre-cleaning	Flush for 5 minutes under cold water, operating moving parts. Clean the instrument with a soft brush (such as MED100.33 from Medisafe GmbH) until no more residues are visible.
Ultrasound and disinfection	Ultrasound bath 35 kHz at room temperature, 10 minutes, cleaning/ disinfection solution 2 % Bomix® plus (Bode Chemie).
Secondary cleaning	Rinse hard to reach areas for 20 seconds with a spray nozzle as needed, then rinse the entire instrument for 30 seconds with demineralized water.

*Proof of fundamental suitability of the instruments for manual cleaning and disinfection was provided for a worst-case product by an independent, accredited test laboratory using the instrument disinfectant Bomix® plus (Bode Chemie) with a concentration of 2 %, test report no. 07015-2 dated 24.11.2015.*

**Machine cleaning and disinfection:**  
 In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

• Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported.  
 • The optionally available storage trays (TAB1-B) ensuring secure support.  
 For improved cleaning results, the lumina of the flushing and suction forceps can be connected to the flushing connection of the CDD using the existing Luer-Lock connections.

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.5 % (5 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 1 minute
Thermal disinfection	90±2 °C, 5 minutes

*Proof of fundamental suitability of the instruments for effective machine cleaning and disinfection was provided by an independent, accredited test laboratory using the Miele G7836 CD disinfecter (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent deconex® 28 ALKA ONE-x with a concentration of 0.5 % (Borer Chemie AG, Zuchwil, Switzerland) (test report no. 111738-10 dated 11.05.2011).*

• Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the instruments and are permitted according to the A<sub>0</sub> concept.

comparable A<sub>0</sub> value>3000. If a different cleaning agent is used, only products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie), for instance regarding the pH value and compatibility with plastics. Contact the responsible supplier or hygiene officer in case of doubt.

**Inspection:**  
 Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

**Maintenance:**  
 None

**Packaging:**  
 Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packaging), or wrap the instrument or tray with the cleaned and disinfected instruments in a cotton cloth and store them together in suitable sterilization containers (TAB1:C, D) meeting the following requirements:  
 • EN ISO/ANSI AAMI ISO 11607  
 • Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)  
 • Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

**Sterilization:**  
 Products must be cleaned and disinfected prior to sterilization.  
 • Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C zzgl. tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	min. 3 minutes
Drying time	30 minutes

*Proof of fundamental suitability of the instruments for effective steam sterilization was provided by an independent, accredited test laboratory, test report no. 111739-10 dated 07.06.2011. Typical conditions in clinics and medical practices as well as the process described above were taken into account.*

△ Do not sterilize in hot air!  
 △ Do not sterilize in STERRAD®!  
 △ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

**Storage / Transportation:**  
 Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging.  
 For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik.

Changes reserved. Current version available at: www.sutter-med.de

STERRAD® is a trademark of Johnson & Johnson, Inc.

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