REF: 360130 – 360219  
REF: 360230 – 360239  
REF: 360701 – 360709  
REF: 368587 – 368687  
REF: 370100 – 379999

0297  
REF: 360226, 360228

Sheet 1/2

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**Product / Use / Disposal:**  
Electrosurgical accessories may only be used and disposed of by trained medical personnel.  
These instructions do not substitute the instructions for the electrosurgical unit and other accessories used.

#### Non-sterile.

⚠ Clean and sterilize before first and subsequent uses..

#### Intended use:

To connect electrosurgical accessories (e. g. electrodes, forceps) to a compatible electrosurgical unit.  
Patient plates to be attached to the patient's body for the deflection of current during electrosurgical applications.

#### Product life:

The life cycle of the product depends on the use and treatment it gets.

#### Prior to use:

⚠ Examine the integrity of the insulation as well as the cleanness and integrity of the instrument.  
⚠ Only use safe and sterilized products. Products must be firmly attached.  
⚠ Only use with compatible products and connectors. Never connect to foreign or unknown current sources.  
⚠ Select the correct interface between the product and the electrosurgical unit/instrument:  
Monopolar – bipolar – patient plate.  
⚠ Connect instrument and cable only to an electrosurgical unit that has been switched off or is in standby mode.  
⚠ Disregarding these instructions may lead to burns and electrical shock.

#### During use:

⚠ Hold plug (not cord) to connect or disconnect cable.  
Do not bend. Do not use rolled-up cables.  
Do not attach cables in an unproper way.  
⚠ Do not use near inflammable or explosive substances.  
⚠ Maximum voltage:  
Monopolar cables: 10000 Vpp  
Bipolar cables: 2000 Vpp  
Monopolar handpieces: 8000Vpp  
**Additional safety information for re-usable patient plates („neutral electrodes“ NE), REF: 360226, 36 0228**

Attach NE **flat to the skin surface** and secure with Velcro. Use two Velcro strips for large NE. Do not place NE over metallic implants. Metallic implants (endoprostheses, coils, etc.) must not be in the current path between the active electrode and NE.  
Do not attach NE on or next to body piercings. Piercings must not be in the current path between the active electrode and NE.  
Do not place NE on scarred tissue, protruding bones or hirsute skin areas. Shave skin area if necessary (but do not place NE over skin that is sore or injured from the shaving). Keep a distance of at least 15 cm between an NE and ECG electrode.  
Do not use electrode gel under the NE. Remove remainders of liquids (e.g. disinfectant) carefully before attaching NE.  
After changing the patient's position, make sure that NE is still fully and correctly attached.  
When using high HF performance (> 100 watts) it is better to use split, e.g. monitored single-use NE.  
Make sure that the long side of NE points to the OR field. See fig. 1 for the position of NE.  
⚠ Disregarding these instructions may pose a risk for the patient (e. g. cause burns).

#### Maintenance and cleaning:

Observe nationally applicable regulations and guidelines.  
Disconnect the cable/handpiece from the instrument/electrosurgical unit.

Do not use aggressive/abrasive cleaners. Immediately, and no later than max. two hours after treatment, clean cables under running water to remove coarse impurities.

#### Machine cleaning and disinfection:

Thermal disinfection for min. 5 minutes at >90°C (>194 °F).  
We recommend cleaners suitable for synthetic materials, e. g. deconex 28 ALKA ONE (Borer Chemie AG, CH).  
Use an approved cleaning program. Follow the manufacturer's instructions for dosage and the program run.  
Use distilled or fully desalinated water for final rinsing cycle.  
Store products safely during machine cleaning/disinfection. Protect from mechanical damage.  
Do not clean together with sharp-edged or pointed objects.  
Ensure that the efficacy of the disinfecter used has been validated (CE marking acc. to DIN EN ISO 15883).

#### Sterilisation:

Only use the following sterilization method: Do not use other sterilization procedures. Only sterilize clean and disinfected products. Use suitable wrapping (DIN EN ISO/ANSI AAMI ISO 11607).  
Use paper/foil wrapping or appropriate sterilization containers.  
Steam-sterilize in the autoclave in accordance with DIN EN 13060 / DIN EN 285.  
Fractioned pre-vacuum sterilization:

Temperature 134 °C (273 °F), 3 minutes; max. Temperature 138 °C (280 °F), max. duration 20 minutes.  
Sterilize cables 36 01 85/86/87 at 134 °C (273 °F) for 20 minutes always.  
Make sure that all items are dry after sterilization!  
⚠ Do not sterilize in hot air.

⚠ Destroy instruments that might have been contaminated with prions (CJD). - Do not reuse!

#### Storage/Transport:

Store in a cool and dry place. Do not expose to direct sunlight. Store and transport in safe containers/packaging.  
Roll cables loosely; do not wind up tightly. Do not bend or fold.

For returns only cleaned and disinfected instruments in sterile packaging will be accepted.

#### Special note:

Cables 37 01 35L/K/S/A/G/R/T/W do not conform to 46.104b of IEC 60601-2-2:2007.  
⚠ Connect these cables to the bipolar outlet of electrosurgical units only. When used properly by trained medical staff, the risk involved is acceptable.

Sutter Medizintechnik cannot be held liable for changes to the product after purchase or deviations from these instructions

Subject to change.  
Updates available under [www.sutter-med.de](http://www.sutter-med.de).

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