Bison™

Bipolare Zange
Bipolar Clamp
Pince bipolaire
Pinzas bipolares
Pinça bipolar
Pinza bipolare
Bipolaire tang
Bipolär Tång
Bipolar fettetang
Bipolaarinen pihti
Kleszczyki bipolarne
Διπολική λαβίδα
Биполярный зажим
جفت تنائي القطب
雙極鉗
양 전극 클램프
兩極鉗子

REF:
73 01 00-XX - 73 01 99-XX
73 05 00
Ø 5 mm:
73 02 00-XX - 73 02 99-XX
73 03 00-XX - 73 03 99-XX
73 04 00-XX - 73 04 99-XX

DEUTSCH	1
ENGLISH	2
FRANÇAIS	3
ESPAÑOL	4
PORTUGUES	5
ITALIANO	6
NEDERLANDS	7
SVENSKA	8
DANSK	9
SUOMI	10
POLSKI	11
Ελληνικά	12
Русский	13
العربية	14
日本語	15
한국어	16
中文	17

TAB A	REF
TAB A1	93 01 00
TAB A2	99 00 50

FIG A	REF
-> FIG A: 1-8	73 01 00-XX – 73 01 99-XX
-> Fig A: 7	73 02 00-XX – 73 02 99-XX
-> Fig A: 2a/2b	73 03 00-XX – 73 03 99-XX
-> Fig A: 8	73 04 00-XX – 73 04 99-XX
-> Fig A: 3-6	73 05 00

Sutter Medizintechnik GmbH
Tullastr. 87
79108 Freiburg/Germany
Tel.: +49 (0) 761 515510
Fax: +49 (0) 761 5155130
E-Mail: info@sutter-med.de
www.sutter-med.de
www.sutter-med.com

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SILVERGlide® is a trademark of Silverglide Surgical Technologies, Inc.

STERRAD® is a trademark of Johnson & Johnson, Inc.

FIG A (Komponenten)

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FIG C (Demontage)

Components (FIG A):

- 1 (Movable) jaws
- 2a Shaft
- 2b Rotary knob with 360 degree rotation
- 3 Push button
- 4 Connector
- 5 Movable part of handle
- 6 Rigid part of handle
- 7 Electrode insert (incl. jaws)
- 8 Insulation tube

Non-sterile. - Clean and sterilize before first and subsequent uses.

Electrosurgical accessories may only be used and disposed of by trained medical staff! These instructions do not substitute the instructions for the electrosurgical unit and other accessories used.

Intended use:

Bipolar clamp for grasping, coagulating and cutting of soft tissue in minimally invasive surgery.

Prior to use:

Assembly (FIG B)

- B1: Slide insulation tube 8 on shaft 2a.
- B2: Push electrode insert 7 gently as far as it will go into shaft 2a.
- B3: Hold rotary knob 2b and secure electrode insert 7 by about six turns as far as it can be turned. The thread must run smoothly and not grate or get caught.

⚠ Do not hold by the movable jaws 1 of the electrode insert 7.

B4: Close jaws; if necessary hold closed with fingertips. Press push button 3 to facilitate the insertion of the shaft. Hold handle as shown in the illustration so that the movable part of the handle 5 will fall into posterior position or can easily be manoeuvred into it.

Leave the movable handle piece in posterior position, do not hold it.

Introduce shaft 2a including the electrode insert 7 straight into the handle. The handle now returns to its home position and push button 3 locks into place.

Rotary knob 2b will now rotate around 360 degrees, and the jaws 1 open and close easily when the handle is squeezed.

Connect the instrument properly only to the bipolar outlet of the electrosurgical unit used.

⚠ Connect the instrument with its cable only to an electrosurgical unit that has been switched off or is in standby mode. Disregarding these instructions may lead to burns and electrical shock.

Examine the integrity of the insulation and the entire instrument as well as the completeness, cleanness and functionality of the product.

- B5: Before using the instrument, perform a short **function test**: Adjust the electrosurgical unit to medium power or 30 – 40 watts, place three layers of wet gauze between the jaws, and activate the unit. When functioning properly, steam will evaporate after 2 – 3 sec.
- ⚠ Do not touch the jaws. They may store heat.
- ⚠ Only use complete, functional, safe and sterilized products!

During use:

- ⚠ Only activate current with tissue grasped between the jaws of the instrument!
- Always use the lowest power setting available to achieve the desired surgical effect.
- Max. voltage: 1000Vp.
- ⚠ The jaws may cause injuries!
 - ⚠ The jaws may store heat after use and cause burns!
 - ⚠ Never place the instrument on or close to the patient!
 - ⚠ Do not use near inflammable or explosive substances!
 - ⚠ The uninsulated parts of the working end of the instrument may carry electricity and cause injuries!
- If needed, clean jaws with sterile solution (Aquadest) and dry. Do not use saline solution!

When using together with irrigation instruments, prefer an irrigation solution that is not electro-conductive.

- ⚠ Undesired activation or movement of laparoscopic electrodes outside the surgeon's field of vision may lead to undesired tissue damage.
- ⚠ Gas embolisms may occur in laparoscopic surgery.
- ⚠ (United States:) Not intended for female sterilization.

Disassembly (FIG C):

- C1: Disconnect the instrument from the cable.
 - C2: Press push button 3 and pull shaft 2a simultaneously out of the handle.
 - C3: Hold rotary knob 2b, take the electrode insert 7 at its front end and unscrew it from shaft 2a.
- ⚠ Do not hold the (movable) jaws 1 of the electrode insert.
- C4: Pull the insulating tube 8 from shaft 2a. The clamp has now been disassembled.

Cleaning and sterilization:

Clean the instrument including the thread of the tip thoroughly immediately after use. Do not allow blood and debris to dry on the instrument.

Manual **Pre**-cleaning:

Use a soft cloth or brush to remove blood and debris.

Use a cleaning brush (TAB A1) to clean the jaws.

Rinse instrument including inner lumen thoroughly with clean tap water after cleaning.

Manual pre-cleaning of the lumen: > 15 sec with water, 3 bars pressure.

=> Do not use aggressive/abrasive cleaners.

Electrode inserts must be cleaned in ultrasound: 40°C (104°F), 15 min. with alkaline cleaner.

Machine cleaning/desinfection:

Use an approved cleaning programme: 70°C (158°F), > 7 min. with cleaner; neutralize 5 min., rinse.

We recommend 0.5% concentrated alkaline cleaner.

No residue may remain.

=> Perform visual control prior to sterilization and verify the integrity of the insulation as well as cleanness and integrity of the entire instrument.

Only sterilize clean/disinfecting instruments. Use paper/foil wrapping or appropriate sterilization containers.

Sterilization in the autoclave (fractioned pre-vacuum):

Temperature 134°C (273°F), 5 min.; max. 135°C (275°F), 18 min.

=> Sterilization at high temperatures and long sterilization times shorten the life cycle of the instrument.

Do not sterilize in hot air.

Do not sterilize in STERRAD®.

=> A certain discolouring of the tips of SILVERGlide® forceps is normal and may easily be removed with a polishing cloth (TAB A2).

For returns only sterilized instruments will be accepted.

⚠ Sutter cannot be held liable for changes to the product after purchase or deviations from these instructions.

Information on applicable standards, regulations and existing validations may be obtained from the manufacturer.

Subject to change.

STERRAD® is a trademark of Johnson & Johnson, Inc.

SILVERGlide® is a trademark of Silverglide Surgical Technologies, Inc.

FIG A (Components)
FIG B (Assembly and function test)
FIG C (Disassembly)

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Änderungen vorbehalten/Alterations reserved · E03 · © Sutter Medizintechnik GmbH