Bison™

Bipolare Zange

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STERRAD® is a trademark of Johnson & Johnson, Inc.

FIG A (Komponenten)

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FIG C (Demontage)

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Components (FIG A):

- 1 Movable jaws
- 2a Shaft
- 2b Rotary wheel with 360 degree rotation
- 3 Pushbutton
- 4 Connector
- 5 Movable part of handle
- 6 Rigid part of handle
- 7 Electrode insert (incl. jaws)
- 8 Insulation tube

Non-sterile. – Clean and sterilize before first and subsequent uses.

Electrosurgical accessories may only be used by trained medical staff. These instructions do not substitute the instructions for the electrosurgical unit and other accessories used.

Intended use:

Bipolar clamp designed for grasping, coagulating and cutting of soft tissue during minimally invasive procedures.

Prior to use:**Assembling the product (FIG B)**

- Slide the insulation tube 8 on shaft 2a.
- Push electrode insert 7 gently as far as it will go into shaft 2a.
- Hold rotary wheel 2b and secure electrode insert 7 by about six turns as far as it can be turned. The thread must run smoothly and not grate or get caught.

⚠ Do not hold by the (movable) jaws 1 of the electrode insert.

-Close jaws: If necessary, hold closed with fingertips. Press push button 3 to facilitate the insertion of the shaft.

Hold handle as shown in the illustration so that the movable part of the handle 5 will fall into posterior position or can be easily maneuvered into it.

Leave the movable part of the handle 5 in posterior position; do not hold it.

Slide shaft 2a, including the electrode insert, straight into the handle. The handle now returns to its home position and the pushbutton locks into place.

The product is correctly mounted when rotary wheel 2b can be turned 360 degrees and the jaws 1 can be opened and closed without difficulty using the handle.

Connect the instrument properly only to the bipolar outlet of the electrosurgical unit used.

⚠ Connect the product only to an electrosurgical unit that has been turned off or is in standby mode. Disregarding these instructions may lead to burns and electrical shock! Examine the integrity of the insulation and the entire instrument as well as its completeness, cleanness, and proper function.

Generally the life cycle of the electrode insert depends on the treatment and applications it gets, and may reach up to max. 20 uses.

Before using the instrument, perform a short function test:

Adjust the electrosurgical unit to medium power or 30 to 40 watts. Place three layers of soft gauze between the jaws and activate the unit. When functioning properly, some steam will evaporate after 2 to 3 secs.

⚠ Do not touch the jaws. They may be hot after activation.

⚠ Only use fully mounted, functional, safe and sterilized products!

During use:

Always use the lowest power setting available to achieve the desired surgical effect.

Max. voltage: 2000 Vpp.

Max. power settings

Erbe VIO:

MODE BIPOLAR CUT, effect 4, max. 100 watts

MODE BIPOLAR SOFT COAG, effect 5, max. 120 watts.

Avoid spark formation around the jaws when working with other generators.

⚠ The jaws may cause injuries.

⚠ The jaws may be hot after use and cause burns!

⚠ Never place instrument on or close to patient.

⚠ Do not use near inflammable or explosive substances.

⚠ The uninsulated parts of the working end of the instrument may carry electric current and cause injuries.

If needed, clean jaws with sterile solution (Aquadest) and dry. Do not use saline solution!

When using together with irrigation instruments, prefer an irrigation solution that is not electroconductive.

⚠ Undesired activation or movement of laparoscopic electrodes outside the surgeon's field of vision may lead to undesired tissue damage.

⚠ Avoid unintended activation of the live metal parts of the jaws during activation (FIG D)

⚠ Gas embolisms may occur in laparoscopic surgery.

Disassembly (FIG C):

C1: Disconnect the instrument from the cable.

C2: Press pushbutton 3 and pull shaft 2a simultaneously out of the handle.

C3: Hold rotary wheel 2b with one hand and the electrode insert 7 at its front end with the other hand. Screw the electrode out of shaft 2a.

⚠ Do not hold the movable jaws of the electrode insert.

C4: Pull the insulation tube 8 from shaft 2a. The instrument has now been disassembled.

Reprocessing

Clean instrument thoroughly and immediately (no later than 60 minutes) after use. Also clean the thread inside the jaws and the inside of the shaft tube.

Do not allow blood and debris to dry on the instrument!

Manual pre-cleaning:

Disassemble instrument.

Use a soft cloth or brush to remove blood and debris.

Use a cleaning brush to clean the jaws.

=> Do not use aggressive / abrasive cleaners.

Rinse instrument thoroughly after cleaning. Manual pre-cleaning of the lumen > 15 secs, 3 bars pressure.

General recommendation for the following cleaning steps:

Only use alkaline or mild alkaline cleaner compatible with synthetics, borer deconex 28 ALKA ONE-x.

The instrument and its parts must be cleaned in an ultrasonic bath: alkaline cleaner, 15 min., 40 °C (104 °F).

Machine cleaning/disinfection:

Choose a washer/disinfector that has been tested in accordance with the current standards (DIN EN ISO 15883).

Procedure:

- 4 min. pre-cleaning in cold water
- 6 min. with detergent (0.5% deconex 28 alka one-x) at 70°C (158 °F)
- 3 min. neutralization in warm water >40°C (104 °F)
- 2 min. rinsing with warm water >40°C (104 °F).

No residue may remain.

We recommend the use of the Sutter storage tray (70 17 40) to avoid damage of the shaft tube and electrode insert in the reprocessing treatment.

=> Perform visual control prior to sterilization and verify the integrity of the insulation as well as the cleanness and integrity of the entire instrument. Disassemble and sterilize instrument in paper/foil wrapping or appropriate sterilization containers.

Sterilization in the autoclave (fractioned pre-vacuum):

Use a steam sterilizer and sterilization procedures that have been tested in accordance with the current standards (DIN EN 13060, DIN EN 285, DIN EN ISO 17665): Temperature 134°C (273 °F), 5 min; max. 135°C (275 °F), 18 min.

Drying 10 min.

=> Sterilization at high temperatures and long sterilization times shorten the life cycle of the instrument.

⚠ Do not sterilize electrode inserts in STERRAD®!

For returns, only sterilized instruments will be accepted.

⚠ Sutter Medizintechnik GmbH cannot be held liable for changes to the product after purchase and for deviations from these instructions.

Information about applicable standards, regulations, and existing validations may be obtained from the manufacturer.

Subject to change.

STERRAD® is a trademark of Johnson & Johnson, Inc.

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11/2014 © Sutter