A B


TO-Bite™

Sheet 1/2

REF:
70 09 60SG

CE 0297

Sutter Medizintechnik GmbH
Tullastr. 87
79108 Freiburg/Germany
Tel.: +49 (0) 761 515510
Fax: +49 (0) 761 5155130
E-Mail: info@sutter-med.de
www.sutter-med.de
www.sutter-med.com

REF 89 90 22
01/2015
© Sutter Medizintechnik GmbH

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Use an aldehyde-free disinfectant agent with approved efficiency and suitable for disinfection of MIS

STERRAD® is a trademark of Johnson & Johnson.

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Durig Use:

Always use the lowest power setting available to achieve the desired surgical effect.

△ Max. voltage: 1000 Vpp.

Frequently clean the tips from blood and debris.

△ Forceps tips can cause injuries.

△ Forceps tips may be hot after use and cause burns.

△ Never place forceps on or close to the patient.

△ Do not use with inflammable or explosive substances!

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• Forceps tips may be hot after use and cause burns.

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- EN ISO/ANSI AAMI ISO 11607

- Suitable for steam sterilization: heat resistant up to min. 141 °C (286 °F) and sufficiently permeable to steam.

- Sufficient protection of the instruments / sterilization packaging from mechanical damage.

Manual pre-treatment, machine cleaning/ disinfection, packaging, sterilization

Observe nationally applicable regulations and guidelines.

Disconnect cable from instrument.

Do not allow blood and debris to dry on the instrument.

Remove with a soft cloth or brush.

Do not use sharp/abrasive cleaners.

△ Do not put in hydrogen peroxide (H₂O₂)!

Always use machine cleaning /disinfection to treat the instrument.

Manual pre-cleaning:

• Clean instrument thoroughly with cold water after use (within max. one hour) to remove all visible impurities. Use a soft brush for the manual removal of impurities; never use aggressive or abrasive cleaners, steel wool or metal brushes. Mind particularly the tips of the instruments as well as the Luer connectors [1] and the inner lumen of the suction channel.

• Thoroughly rinse the channels with a cleaning pistol for min. 10 sec. while closing the suction interruptor [2] with a finger.

• Clean the contaminated channels under running water with an appropriate cleaning brush [3] and make sure that the openings at the end are not congested. The tip of the cleaning brush must come out at the end of the suction channel.

• Thoroughly rinse the channels again with a cleaning pistol for min. 10 sec. while closing the suction interruptor [2] with a finger.

• Proceed to clean the instrument in ultrasound: 40 °C (104 °F), 15 minutes, mild alkaline cleaner, 0,5 % concentration, deconex 28 Alka one-x (Borer Chemie).

=> Sterilization at high temperatures and long sterilization times shorten the life cycle of the instrument.

Ensure that the instrument is dry after sterilization.

△ Do not sterilize in hot air!

△ Do not sterilize in STERRAD®!

△ Destroy instruments that might have been contaminated with prions (CJD). – Do not reuse!

Steam sterilize in the autoclave.

• Fractionated pre-vacuum (ensure that product is sufficiently dry).

• Steam sterilization according to EN 13060 / EN 285

• Validation according to EN 17665 (previously: EN 554/ANSI AAMI ISO 11134)

• Max. temperature 138 °C (280°F) plus tolerances according to EN ISO 17665.

• Sterilization time (exposure time at sterilization temperature) with triple-fractionated pre-vacuum for min. 3.5 minutes at 132 °C (270 °F) or 18 minutes (deactivation of prions).

• Drying time 20 minutes.

The fundamental suitability of the instruments for effective steam sterilization was demonstrated by an independent accredited test laboratory using steam sterilization device Autoclave 6-6-6 Selectomat HP.

=> Sterilization at high temperatures and long sterilization times shorten the life cycle of the instrument.

Ensure that the instrument is dry after sterilization.

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