REF:
704462
704495
704499

CE 0297

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EN

Product / User / Disposal:

Electrosurgery accessories may only be used and disposed of by qualified medical personnel!

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.

Sterile product:⚠ Do **NOT** recondition the instrument.The thermolabile plastics used have **not** been validated regarding potentially suitable cleaning and sterilization processes.

Dispose of the instrument according to internal hospital guidelines for sharp, biologically contaminated items.

Intended use:

The bipolar electrodes are intended for coagulating or cutting soft tissue during surgical procedures. They may only be used with RF generators from the company Sutter Medizintechnik GmbH!

Service life:For single use only.**Before using:**

⚠ Prior to use, inspect the product for:

- Sound condition of the sterile packaging
- Integrity of the instrument

Do **NOT** use the instrument in the following cases:

- Discernible damage to the sterile packaging or electrode
- The package has gotten damp or was opened
- The sterility indicator does not have the required color
- The expiration date has been exceeded
- Always have a sterile bipolar electrode on hand as a replacement!

During use:

⚠ Only connect the RaVoR™ single-use to RF generators (CURIS®, BM-780 II) from Sutter Medizintechnik GmbH when these are shut off or respectively not activated („standby“). Failure to comply can lead to burns or electric shocks for the user or patient.

⚠ The RaVoR™ single-use can be connected directly to the CURIS® 4 MHz radio frequency generator.


⚠ The RaVoR™ single-use are encoded for use with the CURIS® 4 MHz radio frequency generator. After connecting the cable plug to the RF generator, the following program indicator has to light up briefly before the program setting is accepted:

- REF 70 44 62: IC 01 EN
- REF 70 44 95: IC 01 EN
- REF 70 44 99: IC 02 EN

⚠ The RaVoR™ single-use can only be connected to the BM-780 II radio frequency generator using the adapter (REF 37 38 51).

⚠ When using the RaVoR™ single-use with the BM-780 II radio frequency generator, the following settings must be noted:

- Setting on the power adjuster (dial): 2-4
- Application time: 5-15 s depending on indication

⚠ Never lay down the instrument on or in the immediate vicinity of the patient during surgery. Do NOT put any heavy objects on top.

⚠ Directly after turning off the RF current, electrode tips may be so hot they cause burns.

⚠ Not for use in the presence of combustible or explosive substances!

⚠ Unintentional activation or movement of the activated electrode outside the field of vision can lead to unidentified tissue damage.

⚠ Max. voltage 600 Vp

Storage / Transportation:

⚠ Note the sterility expiry date!

Store in a cool, dry place. Protect from sunlight. Store at 5 °C to +35 °C, relative humidity 10-99 %, non-condensing.

Please note:

Serious incidents that occur in relation to the product must be reported to the manufacturer and the responsible public authority in the member state in which the user and/or patient resides.

Sutter Medizintechnik GmbH expressly warns against modifying the instrument. Any change or deviation from the instructions provided here waives the liability of Sutter Medizintechnik GmbH.

Changes reserved. Current version available at www.sutter-med.de.

Explanation of symbols:

	Observe instructions
	Protect against moisture
	Protect against sunlight
	Single patient use – do not reuse
	Sterilized with ethylene oxide. Sterility guaranteed unless package opened or damaged.
	Temperature limitation

	Manufacturer
	Caution: According to federal law in the USA, the sale of this instrument is restricted; it may only be sold by a doctor or with a prescription
	Attention: Observe important safety-related information!
	May be used until
	Do not use if package is damaged
	Sterile barrier system
	Product is not made with natural rubber latex
	Medical Device

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