Monopolar Electrodes **ARROWtip™** Single-Use

REF: 36 44 20 - 36 44 75

Sutter Medizintechnik GmbH
Tullastr. 87
79108 Freiburg/Germany
Tel.: +49 (0) 761 515510
Fax: +49 (0) 761 5155130
E-Mail: info@sutter-med.de
www.sutter-med.de

Only for US
E-Mail: info@sutter-usa.com
www.sutter-usa.com

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Product / User / Disposal:

Caution: Electrosurgery accessories may only be used and disposed of by qualified medical personnel!

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.

Sterile product:

Caution: Do **NOT** reprocess the instrument.

The thermolabile plastics used have not been validated regarding potentially suitable cleaning and sterilization processes.

Dispose of the instrument according to internal hospital guidelines for sharp, biologically contaminated items.

Intended use:

The ARROWtip™ Single-Use monopolar electrode is an electrosurgical instrument for precision soft tissue dissection. It is a single-use device intended for cutting, dissecting and cauterizing soft tissue. The ARROWtip™ Single-Use monopolar electrode is not intended for use in the central nervous system or in the central circulatory system.

Electrosurgery cables and handles:

The monopolar electrodes have a standard 2.4 mm (3/32") diameter connection and are compatible with monopolar handpieces and RF generators common in the market.

Service life:

For single use only.

Before use:

Warning: Prior to use, inspect the product for:

- Sound condition of the sterile packaging
- Integrity of the instrument
- To connect the electrode, the electrosurgery device must be turned off or in standby mode

- Insert the connection end of the electrode fully and properly into the electrosurgery handpiece
- Pull the protective cap off the end of the electrode

Warning: Failure to comply may lead to burns and electric shocks!

Verify that the neutral electrode is applied correctly!

Do **NOT** use the instrument in the following cases:

- Damage to the sterile packaging or electrode
- The package has gotten damp or was opened
- The sterility indicator does not have the required color
- The expiration date has been exceeded

Always have a sterile monopolar electrode on hand as a replacement!

During use:

- Always work with the lowest output setting for the desired surgical effect.
- Never lay down the instrument on or in the immediate vicinity of the patient during surgery. Do NOT put any heavy objects on top. Lay out cables and store unused instruments so they are isolated from the patient.
- **Warning:** Electrode tips are sharp and may cause injuries!
- **Warning:** Directly after turning off the RF current, electrode tips may be so hot they cause burns.
- **Warning:** Unintentional activation or movement of the activated electrode outside the field of vision can lead to unidentified tissue damage.
- **Warning:** Not for use in the presence of combustible or explosive substances!
- When irrigating instruments are used simultaneously with the electrode, use non-conductive irrigation fluid if possible!
- **Caution:** Max. voltage 3000 Vp

Storage / Transportation:

Caution: Note the sterility expiry date!

Store in dry place. Protect from sunlight. Store at 5 °C to 40 °C, relative humidity 10-90 %, non-condensing.

Sutter Medizintechnik GmbH expressly warns against modifying the instrument. Any change or deviation from the instructions provided here waives the liability of Sutter Medizintechnik GmbH. Changes reserved. Current version available at www.sutter-usa.com.

Explanation of symbols:

	Observe instructions
	Protect against moisture
	Protect against sunlight
	Single use only – not for reuse
	Sterilized with ethylene oxide. Sterile as long as the packaging is unopened and intact
	Temperature limit
	Manufacturer
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Caution: Note accompanying documents
	May be used until
	Do not use if packaging is damaged

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