REF: 700960SG

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STERRAD® is a trademark of Johnson & Johnson, Inc.

English

Product / User / Disposal:

Electrosurgery accessories may only be used and disposed of by qualified medical personnel!

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.

△ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

Intended use:

Bipolar coagulation of soft tissue. Depending on the instrument's scope of functionality, also for suctioning fluids during surgical procedures.

Service life:

A minimum of 20 reconditioning cycles can be expected with proper use.

Before using:

△ Inspect the product for cleanliness, mechanical function, and intact insulation before each use.

We recommend checking the insulation with a suitable test device.

△ Only use sterilized products that are in flawless condition!

A certain discoloration of the non-stick instrument tip is normal and harmless.

To connect the instrument and cables, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

During use:

△ Always work with the lowest output setting for the desired surgical effect.

△ Maximum allowable voltage 500 Vp.

△ Regularly wipe blood and tissue residue from the tips.

△ Instrument tips may cause injuries!

△ After application, instrument tips may be so hot they cause burns!

△ Never lay down instruments on or in the immediate vicinity of the patient!

△ Not for use in the presence of combustible or explosive substances!

Reconditioning:

General information:

Observe national guidelines and regulations!

Disconnect the instrument from the cable!

Do not allow blood and tissue residues to dry up!

Remove blood and tissue residues with a soft cloth or brush!

Do not use sharp tools / scouring agents!

△ Do not immerse in hydrogen peroxide (H₂O₂)!

Manual cleaning and disinfection:

△ The instrument must always be reconditioned by machine – no manual cleaning!

According to the recommendations of the DGSV ("Deutsche Gesellschaft für Sterilgutversorgung", German Association for Sterile Goods Supply), To-BITE™ is assigned to risk group B*. Machine cleaning is required for such products.

"This classification was assigned according to the DGSV 2013 flow chart for the classification of medical devices based on the recommendations of the KRINKO ("Krankenhaushygiene und Infektionsprävention Kommission", Hospital Hygiene and Infection Prevention Commission)/BfArM ("Bundesinstitut für Arzneimittel und Medizinprodukte", Federal Institute for Drugs and Medical Devices), Bundesgesundheitsblatt (Federal Health Bulletin) 2012; 55:1244-1310

Manual pre-cleaning:

• Promptly after use (within 1 hour at most), thoroughly clean the instrument with cold water until all visible contaminants have been removed. Use only a soft brush for manually removing contaminants. Never use aggressive or abrasive cleaners, steel wool or metal brushes. Pay special attention to the instrument tip, the Luer connections [1], and the inner lumen of the suction channel.

• Thoroughly flush the suction channel with a spray nozzle for at least 10 seconds. While doing so, hold the suction interrupt [2] closed with a finger.

• Clean the contaminated suction channel with a suitable cleaning brush [3] under running water, making sure that the openings at the ends are unobstructed. The tip of the cleaning brush must come out end of the suction channel.

• Thoroughly flush the suction channel with a spray nozzle for at least 10 seconds. While doing so, hold the suction interrupt [2] closed with a finger.

• Then clean the instrument in an ultrasound bath: 40 °C, 15 minutes, mildly alkaline cleaner with a concentration of 0.5 %, deconex® 28 ALKA ONE-x (Borer Chemie).

Proof of fundamental suitability of the instruments for manual pre-cleaning was provided by an independent, accredited test laboratory, test report no. 0704041002 dated 02.04.2004.

Machine cleaning and disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883). The CDD must have connections for flushing the instruments.

Process:

• Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported.

The lumina of the instruments must be connected to the flushing connection of the CDD using the Luer-Lock connections provided.

Program steps

Pre-rinse with cold tap water

Clean with 0.5 % deconex® 28 ALKA ONE-x at 70 °C

Neutralize with warm tap water (40-45 °C)

Mid-cycle and/or final rinse with warm tap water (40-45 °C)

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