REF 899002 | 2025-10-29 © Sutter Medizintechnik GmbH

Bipolare Elektroden / Bipolar Electrodes

REF:
700462 – 700499

DE Bipolare Elektroden	5	FI Bipolaariset Elektrodit	39
EN Bipolar Electrodes	9	NO Bipolære Elektroder	43
FR Électrodes Bipolaires	12	PL Elektrody Bipolarne	46
ES Electrodos Bipolares	16	BG Биполярни Електроди	50
IT Elettrodi Bipolari	20	RO Electrozi Bipolari	54
NL Bipolaire Elektroden	24	SK Bipolárne Elektródy	58
PT Elétrodos Bipolares	28	SL Bipolarne Elektrode	62
SV Bipolär Blektrod	32	CZ Bipolární Elektrody	65
DK Bipolære Blektroder	36	HU Bipoláris Elektródák	69

LT	Bipoliariniai Elektrodai	73
LV	Bipolārie Elektrodi	76
EL	Διπολικά Ηλεκτρόδια	80
HR	Bipolarne elektrode	84
RU	Биполярные Электроды	87
SR	Bipolarne elektrode	91
MK	Биполарни електроди	95
TR	Bipolar Elektrotlar	99

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Product / User / Disposal:

Electrosurgery accessories may only be used and disposed of by qualified medical personnel! Dispose of the instrument according to internal hospital guidelines for sharp, biologically contaminated items.

Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory. For further information about electrical safety, we recommend IEC TR 61289 or DIN EN 60601-2-2 supplement 1.

⚠ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

Intended purpose:

Bipolar electrodes for the coagulation / volume reduction / cutting of soft tissue during surgical procedures. The products are intended for all indications, all age groups, and any gender. In case of pregnancy, the usual restrictions that are known to trained medical personnel apply.

Service life:

A minimum of 20 reconditioning cycles can be expected with proper use.

Before using:

- ⚠ Inspect the product for cleanliness, intact insulation, and integrity before each use. Note the sensitive instrument tips in particular.
 - ⚠ Only use sterilized products that are in flawless condition!
- All products must be used with a suitable connection cable from Sutter!
- ⚠ To connect electrodes or accessories, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

During use:

- Operating guidelines for these product groups are available as non-binding instructions. They can be requested from the manufacturer.
- ⚠ Always work with the lowest output setting for the desired surgical effect.
 - ⚠ Electrode tips may cause injuries!
 - ⚠ After application, electrode tips may be so hot they cause burns!
 - ⚠ Never lay down the electrode on or in the immediate vicinity of the patient! Lay cables and store unused instruments so they are isolated from the patient.
 - ⚠ Not for use in the presence of combustible or explosive substances!
 - ⚠ Unintentional activation or movement of the activated electrode outside the field of vision can lead to unwanted burns!
 - ⚠ Maximum allowable voltage 500 Vp.

Reconditioning:

Observe national guidelines and regulations!

Disconnect electrodes from accessories / connection cables!

The overall reconditioning process encompasses pre-cleaning, cleaning / disinfection, and sterilization.

⚠ Machine cleaning / disinfection is preferred due to effectiveness and reproducibility!

⚠ Do not immerse in hydrogen peroxide (H₂O₂)!

Pre-cleaning:

Do not allow blood and tissue residues to dry up; rinse thoroughly with cold water after max. 1 hour! Use a soft brush if needed (no wire brush or similar).

Manual cleaning and disinfection:

Cleaning step	Description
Pre-cleaning	Rinse for 5 minutes under cold water and clean with a soft brush (such as MED100.33 from Medisafe GmbH) until no more residues are visible.
Ultrasound and disinfection	Ultrasound bath 35 kHz at room temperature, 10 minutes, cleaning/ disinfection solution 2 % Bomix® plus (Bode Chemie).

Cleaning step	Description
Secondary cleaning	Rinse hard to reach areas for 20 seconds with a spray nozzle as needed, then rinse the entire instrument for 30 seconds with demineralized water.

Machine cleaning and disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).

- Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported.

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.3 % (3 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 2 minutes
Thermal disinfection	90±2 °C, 5 minutes

- Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the instruments and are permitted according to the A_0 concept, comparable A_0 -value > 3000. If a different cleaning agent is used, only use products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie), for instance regarding the pH value and compatibility with plastics. Please contact your supplier or hygiene officer in case of questions.

Inspection:

Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

Maintenance:

None

Packaging:

Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packaging – note sharp instrument tips!), or store the instrument in a suitable sterilization container meeting the following requirements:

- EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
- Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

Sterilization:

Products must be cleaned and disinfected prior to sterilization.

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	min. 3 minutes
Drying time	min. 30 minutes

⚠ Do not sterilize in hot air!

⚠ Do not sterilize in STERRAD®!

⚠ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

Storage / Transportation:

Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging. For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Please note:

The instructions listed above have been validated by the manufacturer as suitable for preparing a medical device for reuse. It is the responsibility of the reconditioner to ensure that actual reconditioning with the equipment, materials, and personnel used in the reconditioning facility achieves the desired result.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The product may only be repaired by the manufacturer or an agent expressly authorized by the manufacturer. Otherwise the warranty and any possible additional liability claims against the manufacturer are voided.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik GmbH.

Changes reserved. Current version available at www.sutter-med.de.

STERRAD® is a trademark of ASP Global Manufacturing GmbH

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