Sutter
Sutter Medizintechnik GmbH
Alfred-Walz-Str. 22
79312 Emmendingen
Telefon: +49 (0) 7641 962 56 0
info@sutter-med.de

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Monopolare HF-Elektroden Art.-Nr.: R360440 (300-0), R360450 (301-0), R360442 (305-2), R360462 (344-2), R360816 (381-0), R360818 (382-0), R360817 (383-0), R360811 (384-0), R360812 (384-1), R360813 (384-2), R360814 (384-3), R360815 (384-4), R360350 (386-1), R360352 (386-2), R360353 (386-3), R360354 (386-4), R360325 (386-7), R360328 (386-8), R360320 (390-0), R360804 (390-1), R360321 (390-5), R360510 (401-0), R360511 (402-0), R360513 (404-0), R360513 (404-0), R360514 (405-0), R360518 (416-0), R360521 (418-0), R360540 (490-0).

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Maximum rated voltage of accessory (U_{max}):

Article No.	U _{max}
R360440 (300-0), R360450 (301-0), R360442 (305-2), R360462 (344-2), R360816 (381-0), R360818 (382-0), R360817 (383-0), R360811 (384-0), R360812 (384-1), R360813 (384-2), R360814 (384-3), R360815 (384-4), R360330 (386-1), R360332 (386-2), R360340 (386-3), R360342 (386-4), R360325 (386-7), R360328 (386-8), R360320 (390-0), R360511 (401-0), R360510 (402-0), R360511 (402-0), R360513 (404-0), R360514 (405-0), R360518 (416-0), R360521 (418-0), R360540 (490-0).	4.3 kV

See also label and catalog information.

When used in combination with other electrical accessories, the maximum rated voltage of application is determined by the accessory with the lowest rated voltage. (Refer to the "Intended Use" section for more details.)

To determine the maximum rated voltage of this product, consult the Instructions for Use (IFU), product label, or the current product catalog.

If you have any doubts or uncertainties, please contact the manufacturer. Before use, carefully read the complete IFU of this product, as well as the IFUs of any accessories and the HF generator being used.

Strict adherence to all requirements, safety notices, and warnings outlined in the respective IFUs is mandatory. This IFU is not intended for users in the USA.

2 Intended Use

The product is intended for open or endoscopic surgery and is used for cutting and coagulating biological tissue.

User:
Only to be used by skilled medical professionals who have been introduced to the product, e.g. using this IFU.

The user is obliged to inform the patient about all possible risks, residual risks and side effects, especially risks related to the application of HF technology.

Indication:
Product intended for use in open or/and endoscopic surgery.

Contraindication:
Product is not intended for direct contact with the heart or the central circulatory or central nervous system.

Application of high-frequency current may interfere with cardiac pacemakers and in vivo heart defibrillators, so affected patients and patients with other electrically conductive implants must consult a cardiologist prior to the intervention.

Connection and Activation:
Prior to application of monopolar products, it has to be safeguarded that an HF neutral electrode has been applied correctly to the patient and that this device is connected properly to the respective HF generator.

The product is connected to the corresponding output of the HF generator via a suitable HF cable or HF handle.

Activation is done using the buttons of the HF handle or with the foot switch of the HF generator.

Combination/Compatibility:
Before use, compatibility of the product with the intended HF handle/HF cable and HF generator has to be verified.

In case of uncertainties, contact the manufacturer of this product or the manufacturer of the used accessory or HF generator.

The frequency of the used HF generator shall not exceed the maximum frequency of 4 MHz and not exceed the maximum voltage of the accessory (see section "Scope").

The monopolar HF electrodes can be operated on generators from the manufacturers listed below:

ERBE	KL Martin	EMED	BOWA
Covidien	ValleyLab	Tekno	Berchold

For the instructions for use, safety information, and warnings in the user manuals or instructions for use of the accessories and HF generator used.

Follow the instructions and warnings in the instructions for use of the HF generator.

It is recommended to use a smoke evacuation system.

IMPORTANT:
Handle with utmost care.

This does not only apply for the duration of the surgery but also for the complete duration of storage, processing and transport as well as during the process of connecting the product with the HF accessory and HF generator.

This applies especially for the thin components and other sensitive areas, e.g. the insulation. Improper use immediately will result in loss of warranty.

Liability for any damages incurred will not be accepted.

3 Safety Notices - WARNING

Make sure that you have read the user manual of the HF generator carefully before using the products.

See these instructions for use and the label for the maximum rated voltage of the product. The product must be cleaned and disinfected as well as sterilized before first use and after each use according to a validated procedure (EN ISO 17665) (see section "(Re)Preparation: Cleaning, disinfection and sterilization").

If there is any uncertainty, contact the manufacturer.

Prior to usage a visual inspection and function test has to be done (see section "Visual Inspection and Function Test").

Before use, visually inspect the product for pressure marks or damage.

In case of damage, deformation or similar is detected on the product, it is not allowed to use the device.

It has to be replaced by a new one.

At least one (1) cleaned, disinfected and sterilized backup product has to be available.

It is the responsibility of the user to determine the appropriate product size and product type according to their professional judgement and based on the patient's specific indication, preferred surgical technique and history, etc.

Prior to use, ensure that the product is firmly inserted in the HF handle/HF cable. This must be done carefully, in order to avoid damages on the product and/or injuries of patient, surgical personnel or third party.

Excessive force can damage the product. Therefore, the product has to be observed during the complete application.

Exclusion:
Do not activate the product as long as it is in contact with metal objects and/or optics. During an endoscopic intervention the patient must not come into contact with grounded metal objects such as e.g. surgical desk frames, instrument trays etc.

Pay attention that no flammable substances are present in immediate vicinity as otherwise there is danger of explosion.

Do not touch the tip of the product while it is activated. After the electrocautery current is switched off, the product tip may still be hot and may cause unintentional burns.

Do not insert or withdraw the tip from the surgical site while the product is activated.

Accidental activation of the product or moving the tip out of the field of view may cause unintended tissue injury.

Do not activate the product continuously for a long period of time. Keep on longer breaks between activation phases.

(Refer to section "General Safety Notices for HF-Technology").

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisioned effect for the respective intervention.

During the surgery, mechanical forces may lead to deformation or signs of wear of the product. If, abrasion of the product is detected during the application, the product has to be exchanged by a new one.

Abrasion, adhesion of tissue, discoloration, sooting etc. do not represent a reason for complaint and do not permit claiming the manufacturer's warranty.

4 General Safety Notices for HF-Technology (Excerpt)

In addition to the acknowledged benefits of HF-Technology, the application includes several risks that have to be attended:

Improper use and non-observance of the IFU can lead to unintentional burns of the patient, user or third party.

Continued further education of the surgical personnel is recommended.

a) Environment
Pay attention that no flammable substances (anesthetics, oxidizing gases, endogenous gases etc.) are present in immediate vicinity as otherwise there is danger of explosion.

Only use non-flammable cleaning or disinfection agents.

All oxygen connections must be tight and leak-proof during the procedure.

b) Patient Positioning and Preparing
Ensure proper patient positioning, i.e. only use insulating surgical drapes that are dry, absorbent and liquid-tight.

Isolate conductive surfaces and contact points towards the patient.

Dry tissues of pulp are necessary to use for skin and breast folds as well as between the extremities.

Prior to application, remove any liquids that potentially accumulated in body cavities.

Only use non-flammable disinfectants.

Do not use alcohol-based tinctures.

Only use non-conductive irrigation fluid, if medicinally possible.

Attend the requirements on irrigation fluid for monopolar and bipolar products.

Prior to application, remove any body jewelry from the patient.

Putting a band aid over the body jewelry is not sufficient!

c) Connections
Prior to application, make sure the product is connected properly to the HF-generator and make sure the correct power and performance setting is adjusted and displayed.

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisioned effect for the respective intervention.

d) HF-Neutral Electrode for Monopolar Application

In case of monopolar application, select a HF-neutral electrode suitable for the patient, apply it correctly and connect it properly with the respective HF-generator.

Follow all instructions for proper application of the HF-neutral electrode and, incl. patient protection and patient monitoring, monitoring of the HF-neutral electrode and all further provisions, safety notices and warnings included in the IFU of the HF-neutral electrode.

e) Patient Reactions
All electrocautery devices potentially can cause muscle stimulation during the application. The design of this product minimizes the risk of this undesirable effect.

Nevertheless, muscle stimulation can lead to an unexpected movement of the patient in the surgical field.

f) Handling HF-Accesory
Make sure the accessory is compatible.

Do not touch the instrument tip during the complete application.

As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface, that is not in contact with the patient.

Never store product on the patient.

Never wrap cables around the patient and never lay cables over the patient.

Unintended activation of the product can lead to unintentional burns or other injuries of patient, user or third party.

Continuous product observation. If wear or other peculiarities are detected, remove product immediately from the patient.

Apply only short activation times.

Keep on longer breaks between activation phases.

Only adjust low power settings.

g) Completeness of the System
At the end of the surgery, confirm completeness of the system.

5 (Re) Processing: Cleaning, Disinfection and Sterilization

5.1. Maximum number of reprocessing cycles
Due to the design, materials used, intended use as well as wear and tear, a maximum limit of permitted reprocessing cycles cannot be determined.

When applied according to the Intended Use, the product undergoes natural wear and tear, considering manner and duration of the application as well as manner and frequency of reprocessing.

Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "Visual Inspection" and "Function Test").

Visual inspection and function test, especially the condition of the complete cable, insulation, port, plug component and pins are decisive for whether the product is allowed to be applied again.

Only undamaged and immaculate products are allowed to be used. In case of uncertainties of the condition of the product or peculiarities, it is not allowed to use the product.

5.2. Time requirements for cleaning and disinfection.
Preparation (for cleaning, pre-cleaning and automated cleaning and disinfection has to be done immediately after the application, however not later than 1 hour after the application.

Avoid idle time.

5.3. After the Application
Clean and disinfect product immediately after the application.

Refer to AKI (German organization for instrument reprocessing).

After application of the product, deposit the product carefully (protection of lifetime of the product).

After application, separate contaminated product and deposit it in a suitable container (deposit means "not dry").

Immediately remove gross stain.

Immediately mark damaged or defective products and sort them out.

Accessory, that do not fit for the sieves of the cleaning and disinfection device (CDD), shall be deposited separately in suitable containers.

Close firmly all disposal-transport containers, in order to avoid drying of stain.

Organize transport of contaminated products the way that contamination of the transport ways and environment also is avoided (closed transport).

Unused reusable products have to be reprocessed as well.

Refer to AKI (German organization for instrument reprocessing), Red Brochure, page 30-32). Note that transport containers must also be cleaned and disinfected after use.

5.4. Validation of (Re)Processing
The following validated processing procedure is recommended.

Equivalent deviant processes are possible.

Then, it is the sole responsibility of the user to safeguard the suitability of the actually applied procedure by suitable means (e.g. validation, routine monitoring, verification of material compatibility etc.).

Automated cleaning and disinfection always are preferable.

The following procedure was validated according to EN ISO 17665.

Additional applicable processing requirements specific to the respective clinical place (operator) as well as national or country specific regulations have to be followed as well.

Never use sharp objects for cleaning.

Disinfectants always have to be rinsed and removed carefully.

5.5. Preparation for Cleaning
Remove product from their packaging.

Place it in a container provided for cleaning.

It is not necessary to disassemble the product.

5.6. Pre-Cleaning
Immediately after the application is completed, pre-clean the product. This however not later than 1 hour after the end of the surgery.

Use tap water (potable water quality) (<40°C) and aldehyde-free, non-fixing disinfectants if applicable.

Thoroughly remove surface stain with a soft brush or synthetic fleece, as otherwise particles or dried secretions may adhere.

This could make subsequent cleaning and sterilization difficult or impossible.

Ensure that areas difficult to reach are cleaned thoroughly and rinsed several times.

Cavities and lumen have to be rinsed intensively using at least 3x20 ml cold tap water (<40°C) with the aids of a rinsing adapter (e.g. from the company Medisafe), with a syringe or with a water jet pistol (>30 Sec.).

This pre-cleaning step always has to be done prior to the manual cleaning or cleaning with the cleaning and disinfection device (CDD).

5.7. Manual Cleaning and Disinfection
Prepare an immersion bath with a suitable fluid cleaning agent.

Use a cleaning agent that is compatible to the disinfectant and suitable for immersion baths.

Follow the instructions of the manufacturer of cleaning agent and disinfectant.

Only use agents suitable for medical devices made from metal and plastics with a pH-value between 5.5 and 12.3.

5.8. Automated Cleaning and Disinfection
Do not use high alkaline cleaning agents. These will impair the lifetime of the product.

• Prepare immersion bath with cleaning agent according to the specific cleaning agent IFU.

• Prepare a separate immersion bath with disinfectant according to the specific disinfectant IFU.

• Immerse the product completely in the ultrasonic bath with cleaning agent (e.g. 0.5% g-igazyme®)

• Clean product in ultrasonic bath using a sonication time of 5 Min. and a frequency of 35 kHz.

• Follow all instructions set forth in the IFU of the cleaning agent, disinfectant and ultrasonic bath.

• Ensure the product will not touch other products or parts in the ultrasonic bath.

• Avoid sonic erosion in the ultrasonic bath are avoided.

• Then, clean product with a soft brush under cold running town water (<40°C).

• Intensively rinse cavities and lumen with a water jet pistol (>30 Sec.) or similar for at least 1 Min.

• Afterwards, rinse product thoroughly for at least 1 min. under running tap water (>40°C) in order to remove any residues of the cleaning agent.

• Inspect product visually on remaining stain.

• In case stain is still present, repeat aforementioned cleaning steps as long as it needs until no soil is present.

• Afterwards, immerse product completely in a disinfectant bath including e.g. Korsolex Plus 3%, for at least 15 min.

• Follow the manufacturer's data for residence time.

• Additionally, Rinse with narrow and areas difficult to access, all cavities and lumen with a syringe several times (at least 3x) using each time 20 ml cold demineralized water.

• Dry product with a lint free wipe and sterile compressed air.

• Dry cavities, lumen and channels with sterile compressed air.

5.9. Automated Cleaning and Disinfection
Only use cleaning and disinfection devices (CDD) with proven efficiency according to EN ISO 15883.

Follow data of the manufacturer of the cleaning and disinfection machine.

Only use agents suitable for medical devices made from metal and plastics with a pH-value between 5.5 and 12.3.

Recommendation: neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG).

Apply program for thermal disinfection.

Follow instructions and data regarding program course and machine.

Product has to be stored safely and protected against mechanical damages during automated cleaning and disinfection.

Do not clean together with sharp edged or pointed objects.

Deposit product in a suitable rinsing basket.

Follow data for loading of the cleaning and disinfection device (CDD).

Use rinsing adapters for products with lumen and connect them according to the instructions in the User Manual of the cleaning and disinfection device (CDD).

Cleaning Program
Start the program course with following parameters:

- 1 min. pre-rinsing with cold water
- Emptying
- 3 min. pre-rinsing with cold water

- Emptying
- 5 min. cleaning at 55°C with 0.5% alkaline cleaning agent
- Emptying
- 3 Min. neutralization with warm tap water (>40°C) and neutralizer (0.1% Neodisher® Z)
- Emptying
- 2 min. interim rinsing with warm demineralized water (>40°C)
- Emptying

Disinfection Program
Automated thermal disinfection considering national requirements regarding AO value (see EN ISO 15883, AO value >3000).

- 5 min. cleaning at 92°C +/-2°C

Drying
30 min. at 90°C

• Remove rinsing adapter

At the end of program course, remove product and inspect it on remaining stain.

In case of residues, repeat automated cleaning and disinfection step as long as it takes until stain is no longer present.

Dry cavities and not sufficiently dried areas with sterile compressed air <2 bar.

Immediately after removal of product and immediately after additional drying on a clean plate, put product in a single use sterilization packing (double packing) from paper or foil or put product in a sterilization container.

Respect requirements for sterilization packaging according to EN ISO 11607 and EN 868.

5.9. Sterilization
Only products that have been cleaned and disinfected are allowed to be sterilized.

Only apply steam sterilization in autoclave (fractioned pre-vacuum with sufficient product drying) for this product.

Adjust sterilization parameters:

- Minimum 134°C and maximum 137°C in saturated steam.
- Holding time at least 5 Min. until max. 20 Min.
- Drying in vacuum for at least 10 Min.
- Sterilizer (Class B) according to valid national standards and regulation (e.g. EN 13060 or EN 285).

Example: Sterilizer Class B, manufacturer: Tuttnauer.

Respect data of the Sterilizer manufacturer regarding load, handling and drying times.

Do not apply hot air, EO-gas, Radiation or Plasma for sterilization, or any other sterilization method for this product.

IMPORTANT:
Prior to usage, let product cool to room temperature.

It is the sole responsibility of the user to maintain the sterile condition of the product after the sterilization process.

In case the aforementioned chemicals and machines for cleaning, disinfection or sterilization are not available, it is the responsibility of the user to validate the procedure actually applied.

Also, if a sterilization method other than described above is applied, this deviating procedure has to be validated by the user accordingly.

5.10. Limitation of Reprocessing
The product life time depending on wear and tear, handling, application time, damages as well as frequency of reprocessing.

Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "Visual Inspection" and "Function Test").

Only an undamaged product is allowed to be reused.

6 Visual Inspection and Function Test

Prior each use, check the entire product, especially the insulation and product tip, for pressure marks and damage.

A product exhibiting damages; pressure points or questionable condition is not allowed to be used and has to be replaced by a new one.

During and after application, tissue may adhere to the product, or sooting may be present on the distal end of the active electrode.

Such adhesions or sooting do not represent a reason for complaint and the product has to be exchanged by a new one.

Due to longer application time, mechanical forces or plasma seam or similar, the product may exhibit deformation or abrasion of the insulation material.

Also, such aspects do not represent a reason for complaint and the product has to be replaced by a new one.

Blockage of the suction channel (if applicable), does not represent a reason for complaint.

Exclusion, that exhibits a blocked suction channel (if applicable), has to be replaced by a new one.

Prior to usage an electrical continuity test has to be done. Furthermore, an insulation test with an insulation test device has to be done.

In case the product does not pass the electrical continuity test and/or insulation test, it is not allowed to use the product any longer and has to be replaced by a new one.