


Featured Product

700952 – Calvian endo-pen®


Qty.	REF	Description
1	700950*	Calvian endo-pen®, 45° angled tip, without suction, working length 180 mm
1	700952*	Calvian endo-pen®, 15° angled tip, without suction, working length 180 mm
1	700954*	Calvian endo-pen®, 70° angled tip, without suction, working length 180 mm
1	700956*	Calvian endo-pen®, forceps tip, 45° tip angled upwards, without suction, working length 180 mm

* All Calvian endo-pen® are available with both EU flat connector or US 2-pin connector.
For more information, please visit: www.sutter-med.de

Forceps jaws (horizontal) 1:1



Duckbill (vertical) 1:1



870010 – CURIS® basic set with single-use patient plates

Qty.	REF	Description	Unit settings / Other accessories
1	360100-01	CURIS® radiofrequency generator (incl. main cord, user's manual and test protocol)	CURIS®
1	360110	Footswitch two pedals for CURIS® (cut & coag), 4 m cable	Calvian endo-pen®: BIPOLAR PRECISE
1	370154L	Bipolar cable for CURIS®, length 3 m	Power adjustment: 15-25 watts
1	360704	Monopolar handpiece (pencil) cut & coag, shaft 2.4 mm, cable 3 m	
1	360238	Cable for single use patient plates, length 3 m	
1 (x50)	360222	Safety patient plates, single use, packing 5 x 10 pcs. (not shown)	



SUTTER MEDIZINTECHNIK GMBH

TULLASTRASSE 87 · 79108 FREIBURG / GERMANY · TEL. +49 (0)761 51551-0 · FAX +49 (0)761 51551-30
WWW.SUTTER-MED.COM · WWW.SUTTER-MED.DE · E-MAIL : INFO@SUTTER-MED.DE