

Calvian endo-pen[®]

Bipolar Forceps



NEW TIP
STYLES



Vertical movement

35°

lateral view

REF 70 09 53
straight-shaped
35° / 1.0 mm tip
working length **18 cm (7.0")**
total length **30 cm (12.0")**

OR
REF 70 09 25
bayonet-shaped
35° / 1.0 mm tip
working length **14 cm (5.5")**
total length **30 cm (12.0")**



Horizontal movement

lateral view

REF 70 09 91
straight-shaped
0.4 mm tip
working length **18 cm (7.0")**
total length **30 cm (12.0")**

OR
REF 70 09 21
bayonet-shaped
0.4 mm tip
working length **14 cm (5.5")**
total length **30 cm (12.0")**

- ✓ Pleasant maneuverability of the endoscope due to
 - slender shaft and handle
 - connector keeps the cable out of the way
- ✓ No assembly, no loose parts
- ✓ Forceps-like handle

- ✓ No retracting tips
 - precise coagulation
 - grasping and manipulating is possible
- ✓ Easy rotation of the instrument
- ✓ Perfect grip due to soft-touch material

Product Information (Ordering)

Each Calvian endo-pen® comes with an instrument tray for safe cleaning and sterilization.

Calvian endo-pen® bipolar forceps

1.0 mm tip width
35° angled tip
straight-shaped
working length: 18 cm (7.0")
total length: 30 cm (12.0")



EU flat connector
70 09 53

with vertical movement



Calvian endo-pen® bipolar forceps

1.0 mm tip width
35° angled tip
bayonet-shaped
working length: 14 cm (5.5")
total length: 30 cm (12.0")



EU flat connector
70 09 25

with vertical movement



Calvian endo-pen® bipolar forceps

0.4 mm tip width
straight-shaped
working length: 18 cm (7.0")
total length: 30 cm (12.0")



EU flat connector
70 09 91

with horizontal movement



Calvian endo-pen® bipolar forceps

0.4 mm tip width
bayonet-shaped
working length: 14 cm (5.5")
total length: 30 cm (12.0")




EU flat connector
70 09 21

with horizontal movement



Instrument tray for one (1) Calvian endo-pen® bipolar forceps (instruments not included)



70 17 78-01

Connection for the rinsing port of the washer-disinfector included

Instrument tray for two (2) Calvian endo-pen® bipolar forceps (instruments not included)



70 17 78-02

Connection for the rinsing port of the washer-disinfector included


Instrument tray for five (5) Calvian endo-pen® bipolar forceps (instruments not included)



70 17 78-05

Connection for the rinsing port of the washer-disinfector included

Instrument tray for ten (10) Calvian endo-pen® bipolar forceps (instruments not included)



70 17 78-10

Connection for the rinsing port of the washer-disinfector included

Product availability is subject to regulatory approval in individual markets. Products may therefore not be available in all markets. The listed lengths and sizes serve as a guideline and may be rounded up or down. The actual lengths may vary slightly.



PRECISION ELECTROSURGERY
Made in Germany

SUTTER MEDIZINTECHNIK GMBH

ALFRED-WALZ-STR. 22 · 79312 EMMENDINGEN/GERMANY
TEL. +49(0)7641-96256-0 · FAX +49(0)7641-96256-30
WWW.SUTTER-MED.COM · INFO@SUTTER-MED.DE