

# Calvian endo-pen<sup>®</sup>


**non-stick** bipolar forceps



Combining the benefits of the Calvian endo-pen<sup>®</sup> bipolar forceps with our established SuperGliss<sup>®</sup> non-stick technology

[Laboratory tests confirm the non-stick properties that last throughout the lifetime of the instrument.\*1]

**New!**



Bipolar coagulation in endonasal and transnasal endoscopic surgery with Calvian endo-pen<sup>®</sup>



lateral view

### non-stick technology

The material specially developed for SuperGliss<sup>®</sup> non-stick bipolar forceps prevents overheating of the tips during coagulation.



“The Calvian endo-pen<sup>®</sup> bipolar forceps are a highly effective means of coagulation and can be safely used for minimally invasive endoscopic endonasal transsphenoidal surgery”

**R. Gerlach, MD**  
Erfurt (Germany)

- ✓ Improved maneuverability of the endoscope due to
  - slender shaft and handle
  - connector keeps the cable out of the way
- ✓ No assembly, no loose parts
- ✓ Forceps-like handle
- ✓ No retracting tips
  - precise coagulation
  - grasping and manipulating is possible
- ✓ Easy rotation of the instrument
- ✓ Perfect grip due to soft-touch material

\*1 Data on file, Sutter Medizintechnik GmbH, Freiburg (Germany)

## Product Information (Ordering)

Each Calvian endo-pen® non-stick comes with an instrument tray for safe cleaning and sterilization.

### Calvian endo-pen® non-stick



0.7 mm tip width  
15° angled tip  
working length: 10 cm

EU flat connector  
**70 09 89 SG**

US pin connector  
**70 09 89 SGS**



### Calvian endo-pen® non-stick



0.7 mm tip width  
15° angled tip  
working length: 18 cm

EU flat connector  
**70 09 88 SG**

US pin connector  
**70 09 88 SGS**




### Instrument tray, holds 1 Calvian endo-pen® (without instruments)

#### 70 17 78-01

For 1 forceps with EU flat connector or  
US pin connector.

Connection for the rinsing port of the  
washer-disinfector included.



Cables for Calvian endo-pen® non-stick bipolar forceps can be found in the Calvian endo-pen® product catalog.

Products shown on this flyer are available in the EU. Availability in other markets may vary.



**SUTTER MEDIZINTECHNIK GMBH**

TULLASTRASSE 87 · 79108 FREIBURG/GERMANY · TEL. +49(0)761-51551-0 · FAX +49(0)761-51551-30

WWW.SUTTER-MED.COM · INFO@SUTTER-MED.DE