



Ordering Information - Bestellinformation

To-BiTE™ non-stick bipolar clamp - non-stick bipolare Zange

- 70 09 60 SG** To-BiTE™ non-stick bipolar clamp
To-BiTE™ non-stick bipolare Zange
- 37 01 54 R** Bipolar (CURIS®) cable for To-BiTE™, length: 3 m
Bipolares (CURIS®) Kabel für To-BiTE™, Kabellänge: 3 m



CURIS® 4 MHz Radiofrequency Generator - 4 MHz Radiofrequenz-Generator

- 36 01 00-01** CURIS® 4 MHz radiofrequency generator (incl. main cord, user's manual and test protocol)
CURIS® 4 MHz Radiofrequenz-Generator (inkl. Netzkabel, Gebrauchsanweisung und Testprotokoll)
- 36 01 10** Foot switch with two pedals for CURIS® (cut & coag) with holding bracket, cable length: 4 m
Zweipedal-Fußschalter CURIS® (cut & coag) mit Haltebügel, Kabellänge: 4 m
- or/oder 36 01 14** Foot switch with two pedals for CURIS® (cut & coag) without holding bracket, cable length: 4 m
Zweipedal-Fußschalter CURIS® (cut & coag) ohne Haltebügel, Kabellänge: 4 m
- 37 01 54 L** Bipolar cable for CURIS®, cable length: 3 m
Bipolarkabel für CURIS®, Kabellänge: 3 m
- 36 07 04** Monopolar handpiece (pencil) cut & coag, shaft 2.4 mm, cable length 3 m
Monopolarer Handgriff, Schaftdurchmesser 2,4 mm, Kabellänge: 3 m
- 36 02 38** Cable for single-use patient plates, cable length: 3 m
Anschlußkabel für Einweg-Neutralelektroden, Kabellänge: 3 m

Available patient plates - verfügbare Neutralelektroden

- 29 00-5** Single-use patient plate, split, for adults and children, PU 20 x 5 pcs.
Einweg-Neutralelektrode, geteilt, für Erwachsene und Kinder, VE: 20 x 5 Stk.
- 36 02 26** Re-usable rubber patient plate
Wiederverwendbare Gummi-Neutralelektrode



Product availability is subject to regulatory approval in individual markets. Products may therefore not be available in all markets.
Lengths for orientation purposes; may vary slightly.

© Sutter Medizintechnik · Subject to change · REF 1010 D - V02 · printed on acid free paper



PRECISION ELECTROSURGERY
Made in Germany

SUTTER MEDIZINTECHNIK GMBH
ALFRED-WALZ-STR. 22 · 79312 EMMENDINGEN/GERMANY
TEL. +49(0)7641-96256-0 · FAX +49(0)7641-96256-30
WWW.SUTTER-MED.COM · INFO@SUTTER-MED.DE

Bipolar Tonsillectomy with To-BiTE™ non-stick bipolar clamp



Bipolare Tonsillektomie mit der To-BiTE™ non-stick bipolaren Zange

To-BiTE™

FOUR FUNCTIONS

„Four functions in one instrument make tonsillectomies safer, faster, easier.“

P. Tolsdorff, MD, Bad Honnef (Germany)

1

Coag
Koagulieren

2

Dissect
Spreizen

3

Aspirate
Saugen

4

Grasp
Greifen



Made in Germany.



Tonsillectomies can be performed very quickly with little, but precise and effective bipolar coagulation while most of the dissection is “cold” or “blunt”. The glossopharyngeal nerve can be identified and spared.¹



Tonsillar Dissection.*
Präparieren der Tonsille.*



Wound immediately after tonsillectomy.*
Wunde unmittelbar nach der Tonsillektomie.*

¹ Tolsdorff P. Radiofrequency Micro-Bipolar Tonsillectomy. ENT & audiology news. 2010; Vol 18 No. 6: 14

*Courtesy of Dr. P.Tolsdorff

*Mit freundlicher Genehmigung von Dr. P.Tolsdorff