



# Swyng<sup>®</sup> non-stick

Bipolar Forceps  
single-use



PRECISION ELECTROSURGERY  
Made in Germany

## non-stick Technology


---

Our proprietary non-stick technology was developed to reduce sticking and charring during coagulation – thanks to optimal heat transfer capabilities.

## Low weight and optimal balance

---

Light-weight aluminum core and innovative product design featuring a dynamic curved shape that follows the natural ergonomics of the hand for improved comfort and ideal balance.



## Two Product Lines

---



### Standard

The pylon guide prevents scissoring of the tips and facilitates the precise coagulation and dissection in your surgical cases. Available with **straight** and **angled tips**.



### Slim


The slim profile was designed for improved visualization of the surgical field, even in narrow operating corridors.

Standard

## Large Product Portfolio

---

- 18 different patterns
- straight and angled tip styles
- three lengths and tip sizes available
- two product lines



---

Product availability is subject to regulatory approval in individual markets. Products may therefore not be available in all markets. The listed lengths and sizes serve as a guideline and may be rounded up or down. The actual lengths may vary slightly.



**Slim**



0.5 mm  
**78 44 105 SL**

**Slim**



1.0 mm  
**78 44 110 SL**

**Slim**



1.5 mm  
**78 44 115 SL**

US 2-pin connector	
Total length	21.0 cm (8.25")
Working length	8.0 cm (3.0")



**Slim**



0.5 mm  
**78 44 405 SL**

**Slim**



1.0 mm  
**78 44 410 SL**

**Slim**



1.5 mm  
**78 44 415 SL**

US 2-pin connector	
Total length	24.0 cm (9.5")
Working length	11.0 cm (4.25")



single-use



**Standard**

0.5 mm  
**78 44 905 ST**

**Standard**

1.0 mm  
**78 44 910 ST**

**Standard**

1.5 mm  
**78 44 915 ST**

US 2-pin connector	
Total length	19.0 cm (7.5")
Working length	6.0 cm (2.5")



**Standard**

1.0 mm  
45° angled up  
**78 44 710 ST**

**Standard**

1.0 mm  
45° angled down  
**78 44 810 ST**

**Standard**

0.5 mm  
**78 44 105 ST**

**Standard**

1.0 mm  
**78 44 110 ST**

**Standard**

1.5 mm  
**78 44 115 ST**

US 2-pin connector	
Total length	21.0 cm (8.25")
Working length	8.0 cm (3.0")



**Standard**

1.0 mm  
45° angled up  
**78 44 610 ST**

**Standard**

0.5 mm  
**78 44 405 ST**

**Standard**

1.0 mm  
**78 44 410 ST**

**Standard**

1.5 mm  
**78 44 415 ST**

US 2-pin connector	
Total length	24.0 cm (9.5")
Working length	11.0 cm (4.25")




**DYNAMIC  
CURVED  
SHAPE**

Our disposable non-stick bipolar forceps line!

## The right instrument for your surgical needs

The Swyng® non-stick bipolar forceps offer precision and versatility in the operating room with a wide range of available models!



© Sutter Medizintechnik · Subject to change · REF 1329 B\_USA – V10 · printed on acid free paper

**SUTTER MEDICAL TECHNOLOGIES USA, INC.**

120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339

WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM